 Orthofix Inc.  
3451 Plano Parkway  
Lewisville, Texas 75056-9453 U.S.A.  
1-214-937-3199  
1-888-298-5700  
www.orthofix.com

Australian Sponsor  
Emergo Australia  
Level 20, Tower II  
Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

 Medical Device Safety Service  
(MDSS)  
Schiffgraben 41  
30175 Hannover  
Germany  
+49 511 6262 8630  
www.mdss.com

---

Device System Name:  
**CETRA<sup>™</sup> Anterior Cervical Plate System**

 Click directory below for desired language

<b>English</b> 	<b>2-3</b>
<b>Español</b> 	<b>4-5</b>
<b>Français</b> 	<b>6-7</b>
<b>Italiano</b> 	<b>8-9</b>
<b>Deutsch</b> 	<b>10-11</b>
<b>Português</b> 	<b>12-13</b>
<b>日本語</b> 	<b>14-15</b>

**Orthofix Inc.**  
 3451 Plano Parkway  
 Lewisville, Texas 75056-9453 U.S.A.  
 1-214-937-3199  
 1-888-298-5700  
 www.orthofix.com

**Australian Sponsor**  
 Emergo Australia  
 Level 20, Tower II  
 Darling Park  
 201 Sussex Street  
 Sydney, NSW 2000  
 Australia

**EC REP** Medical Device Safety Service  
 (MDSS)  
 Schiffgraben 41  
 30175 Hannover  
 Germany  
 +49 511 6262 8630  
 www.mdss.com

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### Description:

The Cetra Anterior Cervical Plate System is comprised of an assortment of non-sterile, single use, titanium alloy (Ti6Al4V ELI per ASTM F136) plates and screws that allow a surgeon to build a temporary anterior cervical implant construct. The plate is attached to the anterior aspect of the vertebral body, by means of screws, to the cervical spine. The system includes the necessary instrumentation to assist in the surgical implantation of the devices.

### Indications for Use:

The Cetra Anterior Cervical Plate System is a temporary implant intended for anterior fixation to the cervical spine from C2 to C7 and indicated for:

1. Degenerative disc disease (defined as back pain of discogenic origin with degenerative disc confirmed by patient history and radiographic studies).
2. Spondylolisthesis.
3. Trauma (i.e., fracture or dislocation).
4. Spinal stenosis.
5. Deformities (i.e., scoliosis, kyphosis, and/or lordosis).
6. Tumor.
7. Pseudoarthrosis.
8. Revision of previous surgery.

### Contraindications:

The Cetra Anterior Cervical Plate System is contraindicated in patients with a systemic infection, with a local inflammation at the bone site, or with rapidly progressive joint disease or bone absorption syndromes such as Paget's disease, osteopenia, osteoporosis, or osteomyelitis. Do not use this system in patients with known or suspected metal allergies. Use of the system is also contraindicated in patients with any other medical, surgical or psychological condition that would preclude potential benefits of internal fixation surgery such as the presence of tumors, congenital abnormalities, elevation of sedimentation rate unexplained by other disease, elevation of white blood cells or a marked shift in white blood cell differential count.

### Potential Adverse Events:

All of the possible adverse events associated with spinal fusion surgery without instrumentation are possible. With instrumentation, a listing of possible adverse events includes, but is not limited to:

1. Early or late loosening of any or all of the components.
2. Disassembly, bending, and/or breakage of any or all of the components.
3. Foreign body (allergic) reaction to implants, debris, corrosion products, graft material including metallosis, straining, tumor formation, and/or auto-immune disease.
4. Pressure on the skin from component parts in patients with inadequate tissue coverage over the implant possibly causing skin penetration, irritation, and/or pain.
5. Post-operative change in spinal curvature, loss of correction, height, and/or reduction.
6. Infection.
7. Vertebral body fracture at, above, or below the level of surgery.
8. Loss of neurological function, including paralysis (complete or incomplete).
9. Non-union, delayed union.
10. Pain, discomfort, or abnormal sensations due to the presence of the device.
11. Hemorrhage.
12. Cessation of any potential growth of the operated portion of the spine.
13. Death.

**Note:** As with any major surgical procedure, there are risks involved in orthopedic surgery. Infrequent operative and postoperative complications known to occur are: early or late infection, which may result in the need for additional surgeries, damage to blood vessels, spinal cord or peripheral nerves, pulmonary emboli, loss of sensory and/or motor function, impotence, permanent pain and/or deformity. Rarely, some complications may be fatal.

### Warnings and Precautions:

1. The Cetra Anterior Cervical Plate System is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical spine.
2. To optimize bony union, perform an anterior microdiscectomy or corpectomy as indicated.
3. To facilitate fusion, a sufficient quantity of autograft or allograft material should be used.
4. Always orient the plate along the midline of the spine.
5. Excessive torque applied to the screws when seating the plate may strip the threads in the bone.
6. Failure to place screws within intended angle range may result in plate locking mechanism fracture or screw back-out past the plate locking mechanism.
7. Failure to achieve arthrodesis will result in eventual loosening and failure of the device construct.
8. Do not reuse implants; discard used, damaged, or otherwise suspect implants.

9. Reuse of the devices labeled as single-use could result in injury or re-operation due to breakage or infection. Do not attempt to re-sterilize single-use implants that come in contact with body fluids.
10. The health care provider must thoroughly clean the instruments prior to steam sterilization.
11. Do not combine implants of dissimilar metals, as it may result in galvanic corrosion.
12. The selection of proper size, shape and design of the implant for each patient is crucial to the safe use of this device in pediatric patients.
13. Potential risks identified with the use of this device system, which may require additional surgery, include: device component fracture, loss of fixation, non-union, and fracture of the vertebra, neurological injury, and vascular or visceral injury.
14. The implantation of anterior cervical plate systems should be performed only by experienced spinal surgeons with specific training in the use of anterior cervical plate systems because this is a technically demanding procedure presenting a risk of serious injury to the patient.
15. Based on fatigue testing results, the physician/surgeon should consider the levels of implantation, patient weight, patient activity level, other patient conditions, etc. which may impact on the performance of the system.
16. The safety, efficacy and performance of the system have been established for conditions in which the system is used as intended and when used as identified in the Indications for Use. Performance of the system has not been evaluated for use that is contrary to the intended use, indications for use or for use that is contraindicated. Failure to use the system as indicated could detrimentally affect the performance of its components.
17. The correct handling of the implant is extremely important. Implants should not be excessively or repeatedly bent, notched or scratched. These operations can produce defects in surface finish and internal stress concentrations, which may become the focal point for eventual failure of the device.

### MRI Compatibility Information:

The Cetra Anterior Cervical Plate System has not been evaluated for safety and compatibility in the Magnetic Resonance (MR) environment. This system has not been tested for heating or migration in the MR environment.

### Cleaning:

Cetra Anterior Cervical Plate System instruments and implants are provided clean but not sterile. Once an implant comes in contact with any human tissue or bodily fluid it should not be re-sterilized or used. Please discard all contaminated implants.

All instruments must be thoroughly cleaned after each use. Cleaning may be done using validated hospital methods or following the validated cleaning processes described below.

None of the instruments in the system require disassembly prior to cleaning.

### From Point of Use:

Whenever possible, do not allow blood, debris or body fluids to dry on instruments. For best results and to prolong the life of the surgical instrument, reprocess immediately after use.

### Preparation for Cleaning:

1. All instruments with moving parts (e.g., knobs, triggers, hinges) should be placed in the open position to allow access of the cleaning fluid to areas which are difficult to clean.
2. Soak the instruments for a minimum of 10 minutes in sterile water prior to the manual or automated cleaning process.
3. Use a soft cloth or a soft plastic bristle brush to remove any visible soil from the instruments prior to manual or automated cleaning. Use a soft plastic bristle brush or a pipe cleaner to remove soil from any inner lumens.

### Manual Cleaning:

1. Upon completing the preparation for the cleaning procedure, prepare the Vesphene® Ise agent per the dilution recommended on the label directions (1 ounce per gallon), or 1 mL of Vesphene® Ise to each 128 mL of potable tap water per the manufacturer's Directions for Use label.
2. Bathe instruments in prepared room temperature solution.
3. Manually agitate instruments in Vesphene® Ise solution for 15 minutes.
4. If visible soil is noted, scrub instruments with a soft plastic bristle brush and use the brush or a pipe cleaner long enough to reach the entire length of any interior lumen(s) to remove the soil.
5. Rinse the instruments in USP <1231> purified water for 1.5 minutes.
6. Hang dry the device.
7. Inspect the instruments for visible soil.
8. If visible soil is noted, repeat the steps listed above.

### Automated Cleaning:

1. Upon completing the preparation for the cleaning procedure, set up the washer/disinfector detergent dose at ½ ounce of Endozime AW Plus® per gallon of water or according to the manufacturer recommendations.

2. If visible soil is noted, scrub instruments with a soft plastic bristle brush and use the brush or a pipe cleaner long enough to reach the entire length of any interior lumen(s) to remove the soil.
3. Place scrubbed instruments into the washer baskets.
4. Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
5. The following automated cleaning cycle is recommended (minimum recommended times are provided for each stage):
  - a. Pre-Wash 1: cold potable water, 2 minutes
  - b. Enzyme/Detergent treatment:
    1. Spray with enzyme/detergent, 20 seconds
    2. Soak, 1 minute
    3. Rinse cold potable water, 15 seconds
    4. Rinse cold potable water, 15 seconds
  - c. Wash  $\geq 65^{\circ}\text{C}$ , 2 minutes using Endozime AW Plus<sup>®</sup>
  - d. Rinse 1: hot potable water, 15 seconds
  - e. Rinse 2: hot potable water, 15 seconds
  - f. Rinse 3: hot potable water, 15 seconds
  - g. Rinse 4: hot potable water, 15 seconds
  - h. Thermal rinse  $\geq 93^{\circ}\text{C}$ , 1 minute
  - i. Heated USP <1231> Purified Water Rinse 1: re-circulating 10 seconds
  - j. Heated USP <1231> Purified Water Rinse 2: non re-circulating 10 seconds
  - k. Dry at  $115^{\circ}\text{C}$ , 7 minutes
6. Inspect the instruments for visible soil.
7. If visible soil is noted, repeat the above listed steps until no visible soil is noted.

**Note:** Certain cleaning solutions such as those containing caustic soda, formalin, glutaraldehyde, bleach, and/or other alkaline cleaners may damage instruments. These solutions should not be used.

**Note:** Visually inspect instruments after cleaning and prior to each use. Discard or return to Orthofix any instruments that are broken, discolored, corroded, have cracked components, pits, gouges, or are otherwise found defective. Do not use defective instruments.

**Sterilization:**

**Sterilization in Orthofix Cases with Blue Wrap:**

The Cetra Anterior Cervical Plate System instruments and implants are supplied NON-STERILE. Prior to use, all instruments and implants should be placed in the appropriate Orthofix case which will be wrapped in a FDA cleared sterilization wrap and placed in the autoclave for sterilization by the hospital using one of the following recommended cycles:

Method: Steam	or:	Method: Steam
Cycle: Gravity		Cycle: Prevac
Temperature: $270^{\circ}\text{F}$ ( $132^{\circ}\text{C}$ )		Temperature: $270^{\circ}\text{F}$ ( $132^{\circ}\text{C}$ )
Exposure time: 15 minutes		Preconditioning: Per manufacturer's settings
Drying time: 30 minutes		Exposure time: 4 minutes
Double wrapped		Drying time: 30 minutes
		Double wrapped

Validation and routine monitoring should be performed per ANSI/AAMI ST79 *Comprehensive guide to steam sterilization and sterility assurance in health care facilities*. Other cycles may be used as long as they comply with the above practices and provide a sterility assurance level of  $10^{-6}$ .

**Packaging:**

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The Cetra Anterior Cervical Plate System instruments and implants are provided in a modular case specifically intended to contain and organize the system's components. The system's instruments are organized into trays within each modular case for easy retrieval during surgery. These trays also provide protection to the system components during shipping. Additionally, individual instruments and implants are provided in sealed poly bags with individual product labels.

**Product Complaints:**

Any Health Care Professional (e.g., customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Orthofix Inc., 3451 Plano Parkway, Lewisville, TX 75056, USA, by telephone at 1-214-937-3199 or 1-888-298-5700 or by e-mail at [complaints@orthofix.com](mailto:complaints@orthofix.com).

**Further Information:**

A recommended operative technique for the use of this system is available upon request from Orthofix at the phone numbers provided above.

**Latex Information:**

The implants, instruments and/or packaging material for the Cetra Anterior Cervical Plate System are not formulated with and do not contain natural rubber. The term "natural rubber" includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

**Caution:** Federal law (USA) restricts these devices to sale by or on the order of a physician.

<b>Rx Only</b>	Federal (U.S.A.) law restricts this device to sale by or on the order of a physician		
	See Instructions for Use		Manufacturer
	Single Use Only Do Not Reuse		Authorized Representative
	Catalogue Number		Serial Number
	Provided Non-Sterile		Lot Number

 Orthofix Inc.  
 3451 Plano Parkway  
 Lewisville, Texas 75056-9453 U.S.A.  
 1-214-937-3199  
 1-888-298-5700  
 www.orthofix.com

Australian Sponsor  
 Emergo Australia  
 Level 20, Tower II  
 Darling Park  
 201 Sussex Street  
 Sydney, NSW 2000  
 Australia

 Medical Device Safety Service  
 (MDSS)  
 Schiffgraben 41  
 30175 Hannover  
 Germany  
 +49 511 6262 8630  
 www.mdss.com

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