 Orthofix Inc.  
3451 Plano Parkway  
Lewisville, Texas 75056-9453 U.S.A.  
1-214-937-3199  
1-888-298-5700  
www.orthofix.com

Australian Sponsor  
Emergo Australia  
Level 20, Tower II  
Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

**EC REP** Medical Device Safety Service  
(MDSS)  
Schiffgraben 41  
30175 Hannover  
Germany  
+49 511 6262 8630  
www.mdss.com

---

Device System Name:  
**SambaScrew® SI Fixation System**

 Click directory below for desired language

<b>English</b> <b>EN</b>	<b>2-3</b>
<b>Español</b> <b>ES</b>	<b>4-5</b>
<b>Français</b> <b>FR</b>	<b>6-7</b>
<b>Italiano</b> <b>IT</b>	<b>8-9</b>
<b>Deutsch</b> <b>DE</b>	<b>10-11</b>
<b>Türkçe</b> <b>TR</b>	<b>12-13</b>
<b>Português</b> <b>PT</b>	<b>14-15</b>
<b>Nederlands</b> <b>NL</b>	<b>16-17</b>
<b>Ελληνικά</b> <b>EL</b>	<b>18-19</b>
<b>日本語</b> <b>JA</b>	<b>20-21</b>
<b>繁體中文</b> <b>ZH-TW</b>	<b>22-23</b>
<b>Tiếng Việt</b> <b>VI</b>	<b>24-25</b>

# INSTRUCTIONS FOR USE

Important Information – Please Read Prior to Use

Rx Only  
CE 0086

Orthofix Inc.  
3451 Plano Parkway  
Lewisville, Texas 75056-9453 U.S.A.  
1-214-937-3199  
1-888-298-5700  
www.orthofix.com

Australian Sponsor  
Emergo Australia  
Level 20, Tower II  
Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

EC REP Medical Device Safety Service  
(MDSS)  
Schiffgraben 41  
30175 Hannover  
Germany  
+49 511 6262 8630  
www.mdss.com

## English EN

### Device System Name:

## SambaScrew® SI Fixation System

### Description:

The SambaScrew SI Fixation System is a temporary, multiple component system consisting of non-sterile instruments and non-sterile, medical-grade titanium (Ti-6Al-4V ELI) implants. The system consists of a 9mm diameter, cannulated screw with multiple orifices on its shaft. The Steinmann Pins, Variable Drill Bit, and Packing Tube are single use devices and should be discarded after use.

### Indications for Use:

The SambaScrew SI Fixation System is intended for fixation of sacroiliac joint disruptions. This fixation device is to only be used in skeletally mature patients.

### Contraindications:

The SambaScrew SI Fixation System is contraindicated for use in patients with:

1. Open wounds, infection, presence of tumor, pregnancy, osteoporosis, certain metabolic disorders affecting osteogenesis, certain inflammatory/neuromuscular conditions, and certain neuromuscular deficits which would place an unusually heavy load on the device during the healing period.
2. The implant is made from Ti-6Al-4V ELI (medical-grade titanium alloy). The fixation implant is contraindicated in any individual with a known or suspected allergy, sensitivity or intolerance to metal.

### Potential Adverse Events:

Potential adverse events include, but are not limited to:

1. Allergic reaction or metal sensitivity to foreign body.
2. Cardiovascular system compromise.
3. Death.
4. Decrease in bone density due to stress shielding.
5. Device bending, disassembly, fracture, loosening, migration and/or retropulsion, or subsidence.
6. Dural tears, neural structure injury.
7. Fixation implant migration with or without bone fracture.
8. Fracture of pelvis or sacrum.
9. Gastrointestinal complications (i.e., ileus or bowel perforation).
10. Hemorrhage.
11. Incisional complications (i.e., dehiscence, hematoma).
12. Infection (Incisional or implant site).
13. Loss of spinal mobility or function.
14. Loosening or fracture of fixation implant.
15. Malfunction of fixation device and/or instruments.
16. Malposition of the fixation device.
17. Migration of fixation implant.
18. Neurological injury/deficit which may range from paresthesias to muscle paralysis, loss of rectal or bladder sphincter control, radiculopathies.
19. Organ, connective tissue or nerve damage.
20. Osteoarthritis.
21. Pain, discomfort or abnormal sensation due to device presence.
22. Persistent low back pain.
23. Reproductive system compromise.
24. Screw back-out or breakage possibly leading to local pain, perforation or irritation of adjacent structures.
25. Sepsis.
26. Urological compromise (i.e., infection/retention).
27. Vascular injury.
28. Wound hematoma.

**Note:** Additional surgery may be necessary to correct some of these anticipated adverse events.

### Warnings and Precautions:

1. The SambaScrew device should only be used by surgeons who have been trained in the use of this device. Information on laboratory and clinical training, as well as additional brochures with a detailed description of proper surgical technique, may be obtained from Orthofix. See the SambaScrew SI Fixation System Surgical Manual for instructions on the implant procedure.
2. Infection may occur immediately following implant fixation or a long time afterwards due to transient bacteremia such as caused by dental treatment(s), endoscopic examination or any other minor surgical procedure. To avoid infection at the implant fixation site, it may be advisable to use antibiotic prophylaxis before and/or after such procedures.

### MRI Compatibility Information:

The SambaScrew SI Fixation System has not been evaluated for safety and compatibility in the Magnetic Resonance (MR) environment. This system has not been tested for heating or migration in the MR environment.

### Cleaning:

Implants are provided clean but not sterile. Once an implant comes in contact with any human tissue or bodily fluid, it should not be re-sterilized or used. Please discard all contaminated implants.

All instruments must be thoroughly cleaned after each use using the manual or automated cleaning procedure listed below.

### From Point of Use:

Whenever possible, do not allow the blood, debris or body fluids to dry on instruments. For best results and to prolong the life of the surgical instrument, reprocess immediately after use.

Cleaning may be done using validated hospital methods or following the validated cleaning process described below.

### Preparation for Cleaning:

1. Keep used instruments moist and do not allow blood and/or bodily fluids to dry on the instruments.
2. Remove gross contaminants with a steady stream of lukewarm/cool water (below 110°F/43°C). Rinse each instrument thoroughly. Do not use saline or chlorinated solutions.
3. Disassemble components to expose all surfaces and clean separately. See following sections for assembly and disassembly instructions.

### Instructions for Disassembly:

1. Parallel Guide Set (PN 11-1109-FD7):  
There are three parts to the Parallel Guide Set: 1. Locking Nut, 2. Sliding Block, and 3. Parallel Guide Frame. The Parallel Guide Set requires disassembly prior to cleaning using the following steps:
  - a. Unthread the Nut from the Sliding Block.
  - b. Remove the Sliding Block from the Parallel Guide Frame.
2. Packing Plunger (11-1109-FD8-1) and Packing Tube (11-1109-FD8-2):  
There are two parts: 1. Packing Plunger and 2. Packing Tube (single use only). The Packing Plunger and Packing Tube require disassembly prior to cleaning using the following steps:
  - a. Remove the Plunger from the Tube.
  - b. Discard the Tube, which is single use only. Do not discard the Plunger.
3. Variable Drill Bit (PN 11-1109-FD9-1) and Adjustable Drill Collar (PN 11-1109-FD9-2):  
There are two parts: 1. Variable Drill Bit (single use only) and 2. Drill Collar. The Variable Drill Bit and Adjustable Drill Collar require disassembly prior to cleaning using the following steps:
  - a. Disengage the Drill Collar from the Variable Drill Bit's shaft.
  - b. Slide the Drill Collar over the Drill Bit Shank to take it off.
  - c. Discard the Variable Drill Bit, which is single use only. Do not discard the Adjustable Drill Collar.

### Automated Cleaning:

Instruments should be cleaned using the following automated cleaning process steps:

TREATMENT	TIME (MM:SS)	TEMPERATURE	CLEANING SOLUTION
Enzymatic Wash	04:00	60°C	Steris® Prolystica® 2X Concentrate Enzyme Presoak (1/4 oz. per fluid gallon)
Wash	02:00	Warm Tap Water	Steris® Prolystica® 2X Concentrate Neutral Detergent (1/4 oz. per fluid gallon)
Rinse	02:00	70°C	N/A
Dry	As Needed	20°C – 80°C	N/A

**Note:** Do not use cleaning solutions containing bleach or formalin since they may damage the instrument.

For ultrasonic cleaning, follow manufacturer's specification for water level, concentration levels of cleaning agent and temperature.

When using a mechanical washer, make sure all instruments stay properly in place and do not touch or overlap each other. Always follow the manufacturer's specifications for automatic washers and use a free-rinsing, low-sudsing detergent with a neutral pH (6.0 – 8.5). Due to variations in water quality, the type of detergent and its concentration may require adjustment for optimal cleaning.

Rinse all instruments thoroughly with tap water, de-ionized or distilled water to remove all traces of debris and cleansing agents. Make sure all internal cavities and ratchets are thoroughly rinsed.

Visually inspect for any remaining soil. If necessary, using a soft brush, gently clean, paying close attention to crevices and other hard to clean areas. Re-clean and rinse as described above, if manual cleaning is performed.

Instruments must be thoroughly dried and all residual moisture must be removed before they are stored. Allow to air dry, or use a soft, lint-free, absorbent towel/cloth to dry external surfaces.

All products should be treated with care. Improper use or handling may lead to damage and possible improper functioning of the device.

**Instructions for Assembly and Inspection:**

1. Parallel Guide Set (PN 11-1109-FD7):  
There are three parts to the Parallel Guide Set: 1. Locking Nut, 2. Sliding Block, and 3. Parallel Guide Frame. The Parallel Guide Set requires assembly after cleaning using the following steps:
  - a. Slide the Sliding Block through the underside of the Parallel Guide Frame with the threads leading.
  - b. Slide the Nut over the Sliding Block's shaft and thread it into place, making sure the Sliding Block seats properly.
2. Packing Plunger (PN 11-1109-FD8-1) and Packing Tube (PN 11-1109-FD8-2):  
There are two parts: 1. Packing Plunger and 2. Packing Tube (single use only). The Packing Plunger and Packing Tube require assembly after cleaning using the following step:
  - a. Insert the Plunger into the Tube.
3. Variable Drill Bit (PN 11-1109-FD9-1) and Adjustable Drill Collar (PN 11-1109-FD9-2):  
There are two parts: 1. Variable Drill Bit (single use only) and 2. Adjustable Drill Collar. The Variable Drill Bit and Adjustable Drill Collar require assembly after cleaning using the following steps:
  - a. Slide the Drill Collar, with the arrow pointing towards the proximal end of the Drill Bit, over the Drill Bit Shank until the desired depth position is engaged.
  - b. Flex the parallel beams to release the Drill Collar from the detent engagement.

After cleaning and before use, inspect the instruments for possible damage, wear or non-functioning parts. Discontinue use and replace all instruments that show wear, cracks, corrosion, pitting, discoloration, etc. Damaged or defective instruments should not be used.

**Sterilization:**

The SambaScrew SI Fixation System implants and instruments are provided NON-STERILE and must be thoroughly cleaned and sterilized prior to each use.

Place the dry instruments in the designated locations in the case (see markings on tray surface that indicate proper storage location in the case).

**Sterilization for Instrument Cases with Blue Wrap:**

It is recommended to double wrap instrument cases in standard central supply FDA cleared wrap prior to steam sterilization using the gravity cycles. A single wrap (or pouch) may be used for the pre-vacuum cycle.

Only sterile products should be placed in the operative field. For a 10<sup>6</sup> sterility assurance level, these products are recommended to be steam sterilized by the hospital using one of the three sets of process parameters below:

**Sterilization Parameters for Use in the United States for Instrument Cases:**

Method: Steam  
Cycle: Pre-Vacuum  
Temperature: 270°F (132°C)  
Exposure time: 4 Minutes  
Drying time: NA – Immediate Use Cycle

The following gravity sterilization cycles are not considered by the United States Food and Drug Administration (US FDA) to be standard sterilization cycles. Users should only use sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization containers) that have been cleared by the US FDA for the following selected gravity sterilization cycle specification.

Method: Steam*	or:	Method: Steam
Cycle: Gravity		Cycle: Gravity
Temperature: 273°F (134°C)		Temperature: 250°F (121°C)
Exposure time: 20 Minutes*		Exposure time: 60 Minutes
Drying time: 15-30 Minutes		Drying time: 15-30 Minutes

**Note:** Because of the many variables involved in sterilization, each medical facility should calibrate and verify the sterilization process (e.g., temperatures and times) used for their equipment.

\*Some non-U.S. Health Care Authorities recommend sterilization according to these parameters in order to minimize the potential risk of transmission of Creutzfeldt-Jakob disease, especially for surgical instruments that could come into contact with the central nervous system. This sterilization cycle has not been validated for prion inactivation.

**Sterilization for Screw Caddy with Blue Wrap:**

The SambaScrew SI Fixation System implants are supplied NON-STERILE. Prior to use, all implants should be placed in the appropriate Orthofix caddy which will be wrapped in a FDA cleared sterilization wrap and placed in the autoclave for sterilization by the hospital using one of the following recommended cycles:

Method: Steam	or:	Method: Steam
Cycle: Gravity		Cycle: Pre-vacuum (minimum 4 pulses)
Temperature: 270°F (132°C)		Temperature: 270°F (132°C)
Exposure time: 15 minutes		Exposure time: 4 minutes
Drying time: 30 minutes		Drying time: 30 minutes

Validation and routine monitoring should be performed per ANSI/AAMI ST79 *Comprehensive guide to steam sterilization and sterility assurance in health care facilities*. Other cycles may be used as long as they comply with the above practices and provide a sterility assurance level of 10<sup>6</sup>.

**Packaging:**

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The SambaScrew SI Fixation System instruments and implants are provided in a modular case and caddy specifically intended to contain and organize the system's components. The system's instruments are organized in a modular case for easy retrieval during surgery. The modular case also provides protection to the system components during shipping. Additionally, individual instruments and implants are provided in sealed poly bags with individual product labels.

**Product Complaints:**

Any Healthcare Professional (e.g., customer or user of this system) who has any complaints or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Orthofix Inc., 3451 Plano Parkway, Lewisville, TX 75056, USA, Telephone: 214-937-3199 or 1-888-298-5700 or via email at [complaints@orthofix.com](mailto:complaints@orthofix.com).

**Further Information:**

A recommended operative technique for the use of this system is available upon request from Orthofix at the phone numbers provided above.

**Latex Information:**

The implants, instruments and/or packaging material for the SambaScrew SI Fixation System are not formulated with and do not contain natural rubber. The term "natural rubber" includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

**Caution:** Federal law (USA) restricts these devices to sale by or on the order of a physician.

<b>Rx Only</b>		Federal (U.S.A.) law restricts this device to sale by or on the order of a physician	
	See Instructions for Use		Manufacturer
	Single Use Only Do Not Reuse		Authorised Representative
	Catalogue Number		Lot Number
	Provided Non-Sterile		

 Orthofix Inc.  
3451 Plano Parkway  
Lewisville, Texas 75056-9453 U.S.A.  
1-214-937-3199  
1-888-298-5700  
www.orthofix.com

Australian Sponsor  
Emergo Australia  
Level 20, Tower II  
Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

 Medical Device Safety Service  
(MDSS)  
Schiffgraben 41  
30175 Hannover  
Germany  
+49 511 6262 8630  
www.mdss.com

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