INSTRUCTIONS FOR USE
Important Information – Please Read Prior to Use

Rx Only
CE 0086

 Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

EC REP Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

Device System Name:
CONSTRUX[®] Mini PEEK Spacer System
Cervical Intervertebral Body Fusion Device
Partial Vertebral Body Replacement (VBR) Device

Click directory below for desired language

English EN	2-3
Español ES	4-5
Italiano IT	6-7
Deutsch DE	8-9
Português PT	10-11
英文 ZH-TW	12-13

INSTRUCTIONS FOR USE

Important Information – Please Read Prior to Use

Rx Only
CE 0086

 Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

 Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

English EN

Device System Name:

CONSTRUX® Mini PEEK Spacer System Cervical Intervertebral Body Fusion Device Partial Vertebral Body Replacement (VBR) Device

Description:

The CONSTRUX Mini PEEK Spacer System is comprised of a variety of implants manufactured from PEEK (Polyetheretherketone), as described by ASTM F2026, with titanium markers as described by ASTM F67. The implants are available in two footprint sizes, a small and a large. The implants are available in various heights in one millimeter increments. The superior and inferior surfaces of the implant have a pattern of ripples to provide increased stability and help prevent anterior/posterior movement of the device.

The CONSTRUX Mini PEEK Spacer System is not intended to be used as a stand-alone device. The CONSTRUX Mini PEEK Spacer System must be used with supplemental fixation. The CONSTRUX Mini PEEK Spacer System is used singly and is implanted using an anterior approach. The CONSTRUX Mini PEEK Spacer System is provided non-sterile and requires sterilization prior to use.

Indications for Use:

When Used as a Cervical Intervertebral Body Fusion System:

The CONSTRUX Mini PEEK Spacer System is indicated for spinal fusion procedures at one or two contiguous levels within the cervical spine (C2-T1) in skeletally mature patients with degenerative disc disease (DDD). DDD is defined as neck pain of discogenic origin with degeneration of the disc confirmed by patient history and radiographic studies.

The CONSTRUX Mini PEEK Spacer System is intended for use with autograft and/or allograft comprised of cancellous and/or corticocancellous bone graft and supplemental fixation (e.g., an anterior cervical plate such as the Orthofix Hallmark Anterior Cervical Plate System).

Patients must have undergone a regimen of at least six weeks of non-operative treatment prior to being treated with the CONSTRUX Mini PEEK Spacer System in the cervical spine.

When Used as a Partial Vertebral Body Replacement (VBR) System:

The CONSTRUX Mini PEEK Spacer System is indicated for use in the thoracolumbar spine (T1-L5) for partial replacement (i.e., partial vertebrectomy) of a diseased vertebral body resected or excised for the treatment of tumors, to achieve anterior decompression of the spinal cord and neural tissues, and to restore the height of a collapsed vertebral body. The CONSTRUX Mini PEEK Spacer System is also indicated for treating fractures of the thoracic and lumbar spine.

The CONSTRUX Mini PEEK Spacer System is designed to restore the biomechanical integrity of the anterior, middle, and posterior spinal column even in the absence of fusion for a prolonged period of time. The Partial VBR Device is intended to be used with autograft or allograft.

The CONSTRUX Mini PEEK Spacer System is intended for use with supplemental fixation. The supplemental fixation system that may be used with the CONSTRUX Mini PEEK Spacer System is the Orthofix Spinal Fixation System (SFS) or the Firebird Spinal Fixation System.

Contraindications:

The CONSTRUX Mini PEEK Spacer System, as with other orthopedic implants, is contraindicated for use in patients with:

1. Active infections in which the use of an implant could preclude adequate and appropriate treatment of the infection.
2. Rapidly progressive joint disease or bone absorption syndromes such as Paget's disease, osteopenia, osteoporosis, or osteomyelitis which may prevent adequate fixation.
3. Conditions that may place excessive stresses on bone and implants, such as severe obesity, pregnancy or degenerative diseases. The decision to use this system in such conditions must be made by the physician taking into account the risks versus the benefits to the patient.
4. Prior fusion at the level to be treated.

Potential Adverse Events:

Potential adverse events include, but are not limited to:

1. Failure of the device to provide adequate mechanical stability.
2. Loss of fixation of the implant.
3. Device component failure.
4. Migration or bending of the device.
5. Loss of bony alignment.
6. Non-union.
7. Fracture of bony structures.
8. Resorption without incorporation of any bone graft utilized.
9. Immunogenic response to the implant materials.

Note: As with any major surgical procedure, there are risks involved in orthopedic surgery. Infrequent operative and postoperative complications known to occur are: early or late infection, which may result in the need for additional surgeries, damage to blood vessels, spinal cord or peripheral nerves, pulmonary emboli, loss of sensory and/or motor function, impotence, permanent pain and/or deformity. Rarely, some complications may be fatal.

Note: When used as a partial VBR device, the CONSTRUX Mini PEEK Spacer System is intended for use in affected vertebral body segments that are equal to or smaller than the size of the device. For larger affected vertebral body segments, a larger device indicated for partial or full VBR is recommended.

Warnings and Precautions:

The surgeon should be aware of the following when using implants:

1. The correct selection of the implant is extremely important. The potential for success is increased by the selection of the proper size, shape and design of the implant. The size and shape of human bones present limiting restrictions on the size and strength of implants. No implant can be expected to withstand the unsupported stresses of full weight bearing.
2. Correct handling of the implant is extremely important. Implants should not be bent, notched or scratched. These operations can produce defects in surface finish and internal stress concentrations, which may become the focal point for eventual failure of the device.
3. Single use only. No surgical implants should be reused. Any implant once used should be discarded. Even though the device appears undamaged, it may already have small defects and internal stress patterns that may lead to fatigue failure.
4. Non-sterile; the CONSTRUX Mini PEEK Spacer System implants and instruments are provided non-sterile, and therefore, must be sterilized before each use.
5. Postoperative care is important. The patient should be instructed in the limitations of the implant and should be cautioned regarding weight bearing and body stress on the device prior to secure bone healing.
6. Reuse of devices labeled as single-use could result in injury or re-operation due to breakage or infection. Do not re-sterilize single-use implants that come in contact with body fluids.

MRI Compatibility Information:

The CONSTRUX Mini PEEK Spacer System has not been evaluated for safety and compatibility in the Magnetic Resonance (MR) environment. This system has not been tested for heating or migration in the MR environment.

Cleaning:

CONSTRUX Mini PEEK Spacer System implants are provided clean but not sterile. Once an implant comes in contact with any human tissue or bodily fluids, it should not be re-sterilized or used. Please discard all contaminated implants.

All instruments must be thoroughly cleaned and sterilized after each use. Reprocess all instruments as soon as it is reasonably practical following use to prevent drying. Cleaning may be done using validated hospital methods or following the validated cleaning process described below.

None of the instruments require disassembly prior to cleaning.

From Point of Use:

Wherever possible, do not allow blood, debris, or body fluids to dry on instruments. For best results and to prolong the life of the surgical instrument, reprocess immediately after use.

Preparation for Cleaning:

1. Any instruments with moving parts (e.g., knobs, triggers and hinges) should be separated or placed in the open position to allow better access of the cleaning fluid to the difficult to clean areas.
2. Soak the instruments for a minimum of 10 minutes in sterile water prior to the manual or automated cleaning process.
3. Use a soft cloth or plastic bristle brush to remove any visible soil from the instruments.

Manual Cleaning:

1. Because none of the instruments can be disassembled, use a soft cloth and/or a soft plastic bristle brush to remove any visible soil from the outside and inside of the instrument.
2. Prepare Vesphene® Ilse agent at the use-dilution recommended by label direction – 1 ounce per gallon or 1 mL of Vesphene® Ilse to 128 mL of potable tap water per manufacturer’s recommendations.
3. Bathe instruments in prepared room temperature solution as recommended by the detergent manufacturer.
4. Manually agitate instruments in Vesphene® solution for 15 minutes.
5. Scrub instruments with a soft plastic bristle brush if visible soil is noted and use a soft plastic bristle brush or pipe cleaner long enough to reach the entire length of the interior lumen to remove the soil.
6. Rinse the instruments in USP <1231> purified water for 1.5 minutes.
7. Hang dry the device.
8. Visually inspect the instruments for visible soil.
9. If visible soil is noted, repeat the steps listed above.

Automated Cleaning:

1. Use a soft cloth or soft plastic bristle brush to remove any visible soil from the outside and inside of the instruments. Use a soft plastic bristle brush or pipe cleaner long enough to reach the entire length of the interior lumen to remove the soil from the inner lumens prior to automated cleaning.
2. Set up the washer / disinfectant detergent dose as ½ ounce of Endozime AW Plus® per gallon of water or according to the manufacturer recommendations.
3. Remove the instruments from the instrument case and place them into the washer baskets.
4. Orient instruments in the automated washer’s carriers as recommended by the washer manufacturer.
5. The following automated cleaning cycle is recommended (minimum recommended times are provided for each stage):
 - a. Pre-Wash 1: cold potable water, 2 minutes
 - b. Enzyme/Detergent treatment:
 1. Spray, 20 seconds
 2. Soak, 1 minute
 3. Rinse cold potable water, 15 seconds
 4. Rinse cold potable water, 15 seconds
 - c. Wash ≥ 65°C, 2 minutes using Endozime AW Plus®
 - d. Rinse 1: hot potable water, 15 seconds
 - e. Rinse 2: hot potable water, 15 seconds
 - f. Rinse 3: hot potable water, 15 seconds
 - g. Rinse 4: hot potable water, 15 seconds
 - h. Thermal rinse ≥ 93°C, 1 minute
 - i. Heated USP <1231> Purified Water Rinse 1: re-circulating 10 seconds
 - j. Heated USP <1231> Purified Water Rinse 2: non re-circulating 10 seconds
 - k. Dry at 115°C, 7 minutes
6. Visually inspect the instruments for visible soil
7. If visible soil is noted, repeat the above listed steps until no visible soil is noted.

Note: Certain cleaning solutions such as those containing caustic soda, formalin, glutaraldehyde, bleach, and/or other alkaline cleaners may damage instruments. These solutions should not be used.

Note: Visually inspect instruments after cleaning and prior to each use. Discard or return to Orthofix if instruments are broken, discolored, corroded, have cracked components, pits, gouges, or otherwise found defective. Do not use defective instruments.

Sterilization:

Sterilization in Orthofix Cases with Blue Wrap:

The CONSTRUX Mini PEEK Spacer System implants and instruments are supplied NON-STERILE. Prior to use, all implants and instruments should be placed in the instrument/implant case which will be wrapped in an FDA cleared sterilization wrap and placed in the autoclave for sterilization by the hospital using one of the following recommended cycles:

Method: Steam	or:	Method: Steam
Cycle: Gravity		Cycle: Prevac
Temperature: 270°F (132°C)		Temperature: 270°F (132°C)
Exposure time: 15 minutes		Preconditioning: 4 pulses
Drying time: 30 minutes		Exposure time: 4 minutes
Double wrapped		Drying time: 30 minutes
		Double wrapped

Sterilization in Rigid Sterilization Containers:

When using rigid sterilization containers, clean, inspect and prepare the rigid sterilization container according to the manufacturer’s instructions.

Select the appropriate rigid sterilization container with either a filtered or solid bottom to properly enclose the Orthofix instrumentation cases (recommended 23¼” long x 11¼” wide container). The following sterilization cycle has been validated:

Method: Steam
Cycle: Prevac
Temperature: 270°F (132°C)
Preconditioning: 4 Pulses
Exposure time: 4 minutes
Drying time: 30 minutes

Note: Rigid sterilization containers with solid bottoms cannot be used in gravity steam cycles.

Validation and routine monitoring should be performed per ANSI/AAMI ST79 *Comprehensive guide to steam sterilization and sterility assurance in health care facilities*. Other cycles may be used as long as they comply with the above practices and provide a sterility assurance level of 10⁻⁶.

Packaging:

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The CONSTRUX Mini PEEK Spacer System instruments and implants are provided in modular cases specifically intended to contain and organize the system’s components. The system’s instruments are organized into trays within each modular case for easy retrieval during surgery. These trays also provide protection to the system components during shipping. Additionally, individual instruments are provided in sealed poly bags with individual product labels.

Product Complaints:

Any Healthcare Professional (e.g., customer or user of this system of products) who has any complaints, or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Orthofix Inc., 3451 Plano Parkway, Lewisville, TX 75056, USA, by telephone at 1-214-937-3199 or 1-888-298-5700 or by email at complaints@orthofix.com.

Further Information:

A recommended operative technique for the use of this system is available upon request from Orthofix at the phone numbers provided above.

Latex Information:

The implants, instruments and/or packaging material for the CONSTRUX Mini PEEK Spacer System are not formulated with and do not contain natural rubber. The term “natural rubber” includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Rx Only	Federal (U.S.A.) law restricts this device to sale by or on the order of a physician	
	See Instructions for Use	REF Catalogue Number
	Orthofix.com/IFU	Manufacturer
	Single Use Only Do Not Reuse	ECREP Authorized Representative
	Provided Non-Sterile	SN Serial Number
LOT	Lot Number	

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