 Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

Australian Sponsor
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

 Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

Device System Name:
NewBridge[®] Laminoplasty Fixation System
Interlaminar Fixation Appliance

 Click directory below for desired language

English 	2
Español 	3-4
Français 	5
Italiano 	6-7
Deutsch 	8-9
Português 	10

Orthofix Inc.
 3451 Plano Parkway
 Lewisville, Texas 75056-9453 U.S.A.
 1-214-937-3199
 1-888-298-5700
 www.orthofix.com

Australian Sponsor
 Emergo Australia
 Level 20, Tower II
 Darling Park
 201 Sussex Street
 Sydney, NSW 2000
 Australia

EC REP Medical Device Safety Service (MDSS)
 Schiffgraben 41
 30175 Hannover
 Germany
 +49 511 6262 8630
 www.mdss.com

English EN

Device System Name:

NewBridge® Laminoplasty Fixation System

Interlaminar Fixation Appliance

Description:

The NewBridge Laminoplasty Fixation System is a device comprised of non-sterile, single use, titanium and titanium alloy components. The specially shaped plates, made of commercially pure (CP) titanium conforming to ASTM F67, are designed to fit the anatomy of a dorsally elevated lamina. The plates have screw holes on both ends, which allow for attachment to the vertebral body, and a screw hole in the center for attachment to the allograft.

The screws, made of titanium alloy (Ti6Al4V ELI, per ASTM F136), are available in a variety of lengths and diameters in order to meet individual anatomical requirements.

The NewBridge Laminoplasty Fixation System must always be used with an allograft.

Indications for Use:

The NewBridge Laminoplasty Fixation System is indicated for use in the lower cervical and upper thoracic spine (C3-T3) after a laminoplasty has been performed. The system holds or buttresses the allograft in place in order to prevent expulsion of the allograft or impingement of the spinal cord.

Contraindications:

The NewBridge Laminoplasty Fixation System is contraindicated:

1. For screw attachment or fixation to the posterior elements of the lumbar spine.
2. For single or two level spondylosis without developmental spinal canal stenosis.
3. Under any direct load bearing conditions.
4. When there is focal anterior compression.
5. When there is isolated radiculopathy.
6. When there is loss of anterior column support resulting from tumor, trauma, or infection.

The NewBridge Laminoplasty Fixation System, as with other metallic orthopedic appliances, is contraindicated for use in patients with active infections in which the use of an implant could preclude adequate and appropriate treatment of the infection. The device is also contraindicated for use in patients with known or suspected metal allergies.

Potential Adverse Events:

Potential adverse events include, but are not limited to:

1. Failure of the device to provide adequate mechanical stability.
2. Loss of fixation of the implant.
3. Device component failure.
4. Migration or bending of the device.
5. Immunogenic response to the implant materials.
6. Nerve damage that may occur as a result of surgical trauma.

Note: As with any major surgical procedure, there are risks involved in orthopedic surgery. Infrequent operative and postoperative complications known to occur are: early or late infection, which may result in the need for additional surgeries, damage to blood vessels, permanent pain and/or deformity.

Warnings and Precautions:

The surgeon should be aware of the following when using metallic implants:

1. The correct selection of the implant is extremely important. The potential for success is increased by the selection of the proper size, shape and design of the implant.
2. The correct handling of the implant is extremely important. Plates should not be bent sharply, reverse bent, notched or scratched. All of these operations can produce defects in surface finish and internal stress concentrations which may become the focal point for potential failure of the device.
3. Single use only. No metallic surgical implant should be reused. Any metal implant once used should be discarded. Even though the device appears undamaged, it may already have small defects and internal stress patterns which may lead to fatigue failure.
4. Non-sterile; the plates, screws and instruments are supplied non-sterile, and therefore, must be sterilized before each use.
5. Postoperative care is important. The patient should be instructed in the limitations of the metallic implant.

6. Reuse of devices labeled as single-use could result in injury or re-operation due to breakage or infection. Do not re-sterilize single-use implants that come in contact with body fluids.
7. The safety and effectiveness of the NewBridge Laminoplasty Fixation System has not been established when implanted in the anterior spinal column.

MRI Compatibility Information:

The NewBridge Laminoplasty Fixation System has not been evaluated for safety and compatibility in the Magnetic Resonance (MR) environment. This system has not been tested for heating or migration in the MR environment.

Cleaning:

All instruments and implants must first be cleaned using established hospital methods before sterilization and introduction into a sterile field. Additionally, all instruments and implants that have been previously taken into a sterile surgical field must first be cleaned using established hospital methods before sterilization and reintroduction into a sterile surgical field. Cleaning can include the use of neutral cleaners followed by a deionized water rinse. All products should be treated with care. Improper use or handling may lead to damage and possible improper functioning of the device.

Sterilization:

The NewBridge Laminoplasty Fixation System components are supplied NON-STERILE. Prior to use, all components should be steam sterilized by the hospital using the following recommended cycle:

Method: Steam	or:	Method: Steam
Cycle: Gravity		Cycle: Prevac
Temperature: 250°F (121°C)		Temperature: 270°F (132°C)
Exposure time: 30 minutes		Exposure time: 8 minutes

Packaging:

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The NewBridge Laminoplasty Fixation System instruments and implants are provided in modular cases specifically intended to contain and organize the system's components. The system's instruments are organized into trays within each modular case for easy retrieval during surgery. These trays also provide protection to the system components during shipping. Additionally, individual instruments and implants are provided in sealed poly bags with individual product labels.

Product Complaints:

Any Healthcare Professional (e.g., customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Orthofix Inc., 3451 Plano Parkway, Lewisville, TX 75056, USA, Telephone: 214-937-3199 or 1-888-298-5700 or via email at complaints@orthofix.com.


Further Information:

A recommended operative technique for the use of this system is available upon request from Orthofix at the phone numbers provided above.

Latex Information:

The implants, instruments and/or packaging material for the NewBridge Laminoplasty Fixation System are not formulated with and do not contain natural rubber. The term "natural rubber" includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

Caution: Federal law (USA) restricts these devices to sale by or on the order of a physician.

Rx Only		Federal (U.S.A.) law restricts this device to sale by or on the order of a physician	
	See Instructions for Use		Manufacturer
	Orthofix.com/IFU		Authorized Representative
	Single Use Only Do Not Reuse		Lot Number
	Catalogue Number		Serial Number
	Provided Non-Sterile		

Orthofix Inc.
 3451 Plano Parkway
 Lewisville, Texas 75056-9453 U.S.A.
 1-214-937-3199
 1-888-298-5700
 www.orthofix.com

Australian Sponsor
 Emergo Australia
 Level 20, Tower II
 Darling Park
 201 Sussex Street
 Sydney, NSW 2000
 Australia

EC REP Medical Device Safety Service
 (MDSS)
 Schiffgraben 41
 30175 Hannover
 Germany
 +49 511 6262 8630
 www.mdss.com

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