ORTHOFIX® INTERNAL FIXATION SYSTEM 2	EN	ORTHOFIX® オーツフィックス 内固定システム 63	JA
SISTEMA DI FISSAZIONE INTERNA ORTHOFIX® 7	IT	ORTHOFIX® 内固定系統 69	ZH
SYSTÈME DE FIXATION INTERNE ORTHOFIX® 12	FR	VNITŘNÍ FIXAČNÍ SYSTÉM ORTHOFIX® 73	CS
ORTHOFIX® SYSTEME INTERNEN FIXATIONS 17	DE	SYSTEM ZESPALANIA WEWNĘTRZNEGO ORTHOFIX® 78	PL
SISTEMA DE FIJACIÓN INTERNA ORTHOFIX® 22	ES	ORTHOFIX® INTERNI FIKSIRNI SISTEM 83	SL
ORTHOFIX® INTERN FIKSATIONSSYSTEM 27	DA	ORTHOFIX® 내고정 시스템 88	KO
SISÄINEN KIINNITY SJÄRJESTELMÄ ORTHOFIX® 32	FI	ORTHOFIX® نم ي ل ادل ا ت ب ث تل ا ز اه ج 93	AR
BESKRIVERSE AV ORTHOFIX® INNVEDIG FSTESYSTEM 37	NO	ORTHOFIX® СИСТЕМА ЗА ВЪТРЕШНА ФИКАЦИЯ 98	BG
ORTHOFIX® SYSTEMEN VOOR INTERNE FIXATIE 42	NL	„ORTHOFIX®“ VIDINĖ FIKSAVIMO SISTEMA 103	LT
SISTEMAS DE FIXAÇÃO INTERNA ORTHOFIX® 47	PT	SISTEM FIKSASI DALAMAN ORTHOFIX® 108	MS
ORTHOFIX® INTERNA FIXERINGSSYSTEM 52	SV	SISTEMUL DE FIXARE INTERNĂ ORTHOFIX® 113	RO
ΣΥΣΤΗΜΑ ΕΣΩΤΕΡΙΚΗΣ ΣΤΑΘΕΡΟΠΟΙΗΣΗΣ ORTHOFIX® 57	EL	ORTHOFIX® İNTERNAL FİKSAŞYON SİSTEMİ 118	TR

Orthofix Srl

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Telephone +39 045 671 9000 - Fax +39 045 671 9380

PQ INF H 09/18

## Instruction leaflet - Please read prior to use

# ORTHOFIX® INTERNAL FIXATION SYSTEM

### DESCRIPTION & INDICATIONS FOR USE

The Orthofix Internal Fixation Systems consist of Intramedullary Nailing Implants for femur and tibia, and Fragment Fixation System Implants. These implants are intended as a means of bone stabilization in the management of fractures and in reconstructive surgery. They are not intended to replace normal healthy bone or to withstand the stresses of full weightbearing, particularly in unstable fractures or in the presence of non union, delayed union or incomplete healing. The use of external supports (e.g. walking aids) is recommended as a part of the treatment.

Intramedullary Nailing Implants are indicated in:

1. Traumatic fractures of the diaphysis provided that the epiphyses are closed, so that two locking screws can be inserted outside the fracture both proximally and distally.
2. Pathological fractures.
3. Re-fractures.
4. Non-unions.
5. Reconstructive surgery.

For 2, 3, 4 and 5 above, the limitations are the same as in item 1.

Fragment Fixation System Implants are indicated in:

1. Fractures
2. Bony ligament avulsions
3. Osteotomies

All Orthofix implants are intended for professional use only. Surgeons who supervise the use of Orthofix implants must have full awareness of orthopaedic fixation procedures as well as adequate understanding of the philosophy of the Orthofix System. To promote the proper use of its fixation system, and establish an effective promotional and training tool, Orthofix has developed several manuals of the relevant information (i. e. general philosophy, surgical application, etc.), called "Operative Techniques". These manuals are available as a complimentary service for surgeons who have adopted the Orthofix system.

If you wish to receive a personal copy of one or more of these manuals, which are available in several languages, please contact Orthofix or its local authorized representative, with a description of the medical device to be used.

### CONTRAINDICATIONS

The Orthofix Internal Fixation System is not designed or sold for any use except as indicated.

Use of the Intramedullary Nailing Implants is contraindicated in the following situations:

1. Active infection.
2. General medical conditions including: impaired blood supply, pulmonary insufficiency (i.e. ARDS, fat embolism) or latent infection.
3. Patients with psychological or neurological conditions, who are unwilling or incapable of following instructions for post-operative care.
4. Foreign body sensitivity. Where material sensitivity is suspected, tests should be made prior to implantation.

Use of the Fragment Fixation System is contraindicated in the following situation:

Foreign body sensitivity. Where material sensitivity is suspected, tests should be made prior to implantation. The presence of active infection.

### WARNINGS & PRECAUTIONS

- 1) Correct selection of the implant model and size is extremely important.
- 2) **The Orthofix Tibial Nailing System:** tibial nails are either solid (8mm and 9mm diameter) or cannulated (10 mm, 11mm, 12mm and 13mm diameter). Smaller diameter nails are inserted without or with minimal reaming. Unreamed nails are strong enough to support a stable tibial fracture in most cases, and

should cause less damage to the bone vasculature. An unreamed nail is therefore recommended for all fractures where the external blood supply of the tibia has been disturbed. These include most open fractures, and closed fractures with soft tissue damage of Tscherne types CII and CIII. The usual diameter for an unreamed nail is 9mm, but an 8mm nail may be needed for smaller diameter bones. Sizes above 9mm will nearly always require some reaming. The position and stability of the fracture is also relevant: an unstable or severely comminuted fracture, or a fracture in the metaphyseal area, may need a larger nail for adequate stabilization, and would therefore need to be inserted after reaming. Reaming and nail insertion should not be performed in the presence of an inflated tourniquet, since this may lead to bone or muscle necrosis, and/or compartment syndrome.

- 3) **Calcaneo-Tibial Nails:** in particular selected cases a retrograde nail may be inserted through the heel and across the ankle joint into the distal tibia, for the purpose of arthrodesing the ankle joint. As with the antegrade nail, adequate reaming should be performed for the diameter of nail chosen, the correct locking screws should be inserted, and the operative technique followed carefully. Adequate countersinking of the locking screw heads should be performed to prevent skin irritation. The nail should be dynamised at the correct time to encourage final union.
- 4) **The Orthofix Femoral Nailing System:** to reduce the driving force on the nail and the potential risk of nail damage or fracture comminution, reaming to a width 2 mm greater than the diameter of the proposed nail is recommended. In young patients it is advisable to widen the proximal 10 cm of the femur only by a further 1-2 mm to accommodate the larger proximal 90 mm of the nail. Thus, for a 12/10 mm nail, the proximal 10 cm should be reamed to 13 mm. This is not usually necessary in older patients with softer bone. The distal locking screws of the nail should not be inserted close to the fracture line, as effective screw fixation may be compromised, causing implant failure or loss of fixation.
  - i) After reaming, but prior to nail implantation, the relevant internal nail template should be used. This enables the surgeon to check the amount of reaming to permit nail insertion without hammering; assess nail length accurately, and identify correct positions for the distal locking screws.
  - ii) In very heavy or uncooperative patients, the 13/11 mm or 13/12 mm nail should be used if adequate reaming can be achieved. Again, a stronger nail may be needed if there is severe comminution of the fracture, with little or no bony integrity.
- 5) **The Retrograde Femoral Nailing System:** this system has been developed for the established indication of inserting a femoral nail from the distal end of the bone, for the fixation of fractures or in selected reconstruction procedures. With the exception of the entry point, the same proviso's are indicated as for the antegrade femoral nailing system. Careful reading of the operative technique is essential prior to carrying out any procedure. It is essential that the distal end of the nail is recessed beneath the bone surface, to avoid subsequent damage to the knee joint. The use of a retrograde nailing system in the femur is not recommended in open fractures, or if there is an open wound in the knee. The inevitable risk of infection is around 5%, and this may lead to a septic arthritis of the knee. Alternative methods of treating the fracture should be considered.
- 6) When using compression locking screws with a retrograde femoral nail to treat intracondylar fractures, particular care should be taken to achieve and maintain the reduction before nail insertion, and to choose the correct locking screw length. The screw should then be locked with a compression nut, taking care not to overtighten.
- 7) When indicated where the soft tissues are thin, the heads of nail locking screws should be recessed with the appropriate tool to prevent skin irritation.
- 8) **For both Femoral and Tibial nails,** a mechanical targeting system is provided as part of the instrumentation. It is very important that the anterior hole for the stabilising rod is drilled in the centre of the bone, over the nail. If any doubt exists about this, the position should be checked with X-rays before drilling the hole. No guarantee is given that this system will work every time.
- 9) In very distal fractures, it is important to ensure that the more proximal of the two or three distal locking screws is at least 1 cm distal to the fracture line.
- 10) It is also important that the proximal locking screws are not inserted close to the fracture line, as effective screw fixation may be compromised, causing implant or fixation failure.
- 11) Fracture distraction for any time should be avoided during the operation, and fractures should never be locked in distraction, since this may cause delayed union and/or compartment syndrome.
- 12) Correct handling of the instrumentation and implant is extremely important. Surgical theatre staff should avoid any notching or scratching of the implants, as these factors may produce internal stresses which may become a focus for possible implant breakage.
- 13) If any of the aluminium instruments (the black parts) are hammered, they will be damaged and must be replaced. Hammering of the nails may only be performed through the locking rods, which are made of hardened steel, or with specifically designed insertion tools.
- 14) Bending of the Fragment Fixation System implant during insertion must be avoided, since this may lead to breakage of the implant. When the chamfer is close to the cortex, the speed of insertion must be reduced.
- 15) Wherever a Kirschner wire or Guide Wire is used to guide a cannulated reamer, drill bit or screw into position, the following precautions should be taken:
  - A) The Kirschner or Guide Wire should always be NEW and should NEVER be reused
  - B) The wire should be checked before insertion to exclude any scratches or bends in the wire; if any are found it should be discarded.
  - C) During the introduction of any instrument or implant over a wire, the surgeon should screen the wire tip as continuously as possible to exclude inadvertently driving the wire further than intended.

- D) During each pass of the instrument or implant, the surgeon should check that there is no bony or other debris built up on the wire or inside the instrument or implant which might cause it to bind on the wire and push it forward.
- 16) It is impossible to clean the inside of a cannulated drill bit adequately to exclude organic or other debris remaining after use. CANNULATED DRILL BITS SHOULD THEREFORE NEVER BE REUSED. THEY ARE SINGLE PATIENT USE ONLY. If they are to be used for a second time on the same patient, the surgeon must check that the drill bit is free from obstruction, by removing it from the power unit and passing a wire through it.
- 17) Even when they are new, we recommend that a wire is passed down a cannulated drill bit prior to use, to check that the lumen is free from obstruction.
- 18) Any fixation device may break if subjected to the increased loading caused by delayed or non-union.
- 19) Careful monitoring of the progress of healing must be undertaken in all patients. If callus is slow to develop, other measures may be required to promote its formation, such as dynamisation of the implant, a bone graft, or exchanging the implant for a larger one. It is also important to avoid continued stabilization with a locked small diameter nail of a fracture with delayed union, after 12 weeks in the tibia or 16 weeks in the femur, because of the risk of fatigue failure of the implant. If no callus is visible on X-ray at this time, a second intervention should be considered.
- 20) Advice on weightbearing after insertion of a retrograde femoral nail.
- A) In intra-articular fractures, there should be complete non weightbearing for the first six weeks.
- B) General advice: partial weightbearing (15 kg) after wound healing, increasing slowly up to full weightbearing when there is radiological evidence of callus formation.
- C) In unstable fractures, toe touching initially, increasing to a maximum of 30 kg, progressing to full weightbearing only when there is radiological evidence of callus formation.
- D) With 10 mm nails, partial (15 Kg) up to maximum of 50% body weight, increasing to full weightbearing only when there is radiological evidence of bridging callus formation.
- E) If possible avoid the situation of an unstable fracture fixed with a 10 mm diameter nail. Where this cannot be avoided, the surgeon should take particular care to ensure that the above recommendations for weightbearing are carried out.
- 21) Implants, nail end caps, locking screws and Fragment Fixation System implants must never be reused.
- 22) Patients should be instructed to report any adverse or unanticipated effects to the treating surgeon.
- 23) Nuclear Magnetic Resonance imaging should not be used in any segment to which an implant is applied.
- 24) Removal of the device: the surgeon should make the final decision whether an implant can be removed. Implant removal should be followed by adequate post-operative management to avoid refracture. If the patient is older with a low activity level, the surgeon may decide not to remove the implant, thus eliminating the risks involved in a second surgical procedure.
- 25) To lock the extractor on to the Fragment Fixation System implant, the knob should be turned counter-clockwise.
- 26) All equipment should be carefully examined prior to use to assure proper working condition. If a component or instrument is believed to be faulty, damaged or suspect, it should NOT BE USED.
- 27) Additional equipment may be required for fixation application and removal such as soft tissue retractors, flexible reaming set, cannulated drills, etc.

## POSSIBLE ADVERSE EFFECTS

1. Non-union or delayed union, which may lead to implant breakage.
2. Metal sensitivity, or allergic reaction to a foreign body.
3. Pain, discomfort, or abnormal sensations due to the presence of the device.
4. Nerve damage due to surgical trauma.
5. Bone necrosis.

**Warning: This device is not approved for screw attachment or fixation to the posterior elements (pedicles) of the cervical, thoracic or lumbar spine.**

## Important

A successful result is not achieved in every surgical case. Additional complications may develop at any time due to improper use, medical reasons or device failure which require surgical reintervention to remove or replace the internal fixation device.

Preoperative and operative procedures including knowledge of surgical techniques and proper selection and placement of the internal fixation devices are important considerations in the successful utilization of Orthofix internal fixation devices by the surgeon.

Proper patient selection and the patient's ability to comply with physician instructions and follow prescribed treatment regimen will greatly affect the results. It is important to screen patients and select optimal therapy given physical and/or mental activity requirements and/or limitations. If a surgical candidate exhibits any contraindication or is predisposed to any contraindication, DO NOT USE Orthofix Internal Fixation implants.

## Materials

The Orthofix Internal Fixation System is comprised of stainless steel and aluminium alloy components. Those components which come into contact with patients are made of surgical grade stainless steel.

## STERILE & NON-STERILE Product

Orthofix provides certain internal fixation devices STERILE while others are provided NON-STERILE. Please review the product label to determine the sterility of each device.

## Sterile

Devices or kits provided STERILE are labelled as follows: **Contents of package are STERILE unless package is opened or damaged. Do not use if package is opened or damaged.**

## Non-Sterile

Unless otherwise noted, Orthofix internal fixation components are provided NON-STERILE. Orthofix recommends that all NON-STERILE components be properly cleaned and sterilized following the recommended cleaning and sterilization procedures.

Product integrity and performance are assured only if the packaging is undamaged.

## Cleaning & Maintenance (Attention, see instructions PQ ISP)

Prior to use, NON-STERILE product must be cleaned using a mixture of 70% medical grade alcohol and 30% distilled water. After cleaning, the device and/or system components should be thoroughly rinsed in sterile distilled water and dried using clean non-woven fabric.

To prevent corrosion, the components must be kept dry and detergents with fluoride, chloride, bromide, iodide or hydroxyl ions must be avoided when cleaning, as they will damage the black anodised coating on any Orthofix products, and this may initiate the process of stress corrosion.

Before sterilization, all components should be inspected, since damage to the surface of metal components can reduce strength and fatigue resistance, and may lead to corrosion. If components are damaged in any way, they should be exchanged immediately for new ones.

Assembly of the intramedullary fixation system should then be carried out to check all components.

## Sterilization

The recommended, validated sterilization cycle is:

Methods	Cycle	Temperature	Exposure Time
Steam	Pre-Vacuum	132° - 135°C [270° - 275°F]	Minimum 10 minutes

The instrumentation can be sterilized in the assembled state as long as locking screws are left untightened.

Sterility cannot be assured if the sterilization tray is overloaded. Do not overload the sterilization tray or include additional implants or instruments from any source.

**Orthofix is only responsible for safety and effectiveness for the initial patient use of the internal fixation devices. The institution or practitioner bears full responsibility for any subsequent uses of the devices.**

**CAUTION:** Federal Law (USA) restricts this device to sale by or on the order of a physician.

## RISKS DUE TO THE RE-USE OF “SINGLE USE” DEVICE

### **IMPLANTABLE DEVICE\***

The “SINGLE USE” implantable device\* of Orthofix is identified through symbol  reported on the product label.

After the removal from the patient, the implantable device\* has to be dismantled.

The re-use of implantable device\* introduces contamination risks for users and patients.

The re-use of implantable device\* can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

(\*): Implantable device

Any device which is intended:

Any device intended to be totally / partially introduced into the human body through surgical intervention and intended to remain in place after the procedure for at least 30 days is also considered an implantable device.

### **NON IMPLANTABLE DEVICE**

The “SINGLE USE” non implantable device of Orthofix is identified through symbol  reported on the label or are indicated in the “Instructions For Use” supplied with the products.

The re-use of “SINGLE USE” non implantable device can not guarantee the original mechanical and functional performances, compromising the effectiveness of the products and introducing health risks for the patients.

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**EN** All Orthofix internal and external fixation products should be used together with their corresponding Orthofix implants, components and accessories. Their application should be performed with the specific Orthofix instrumentation, following carefully the surgical technique recommended by the manufacturer in the appropriate Operative Technique Manual.

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