EN	Evolution-TL and Evolution-C Removal Kits - Instructions for Use	SK	Otočné TL a otočné C extrakčné súpravy Návod na použitie
FR	Trousse d'extraction Evolution-TL et Evolution-C - Mode d'emploi	SV	Evolution-TL och Evolution-C borttagningssatser - Bruksanvisning
DE	Evolution-TL und Evolution-C Entnahmekits - Gebrauchsanweisung	IS	Evolution-TL og Evolution-C flutningsúnaður - Notkunarlæðingarn
ES	Kits de extracción Evolution-TL y Evolution-C - Instrucciones de uso	BG	Комплекти за разглобяване Evolution-TL и Evolution-C Инструкции за употреба
IT	Kit di rimozione Evolution-TL ed Evolution-C - Istruzioni per l'uso	HR	Evolution-TL i Evolution-C kompleti za odstranjenje - Upute za korištenje
PT	Kits de Remoção Evolution-TL e Evolution-C - Instruções de utilização	SR	Evolution-TL и Evolution-C комплекти за odstranjivanje Упутства за коришћење
CS	Soupravy Evolution-TL a Evolution-C Removal Kits - Návod k použití	SL	Evolution-TL in Evolution-C kompleta za odstranjevanje - Navodila za uporabo
DA	Evolution-TL og Evolution-C Removal Kits - (Fjernelseskits) Brugsanvisning	RO	Etuse de extracție Evolution-TL și Evolution-C - Instrucțiuni de utilizare
NL	Evolution-TL- en Evolution-C-verwijderkits - Gebruiksaanwijzing	RU	Наборы экстракторов Evolution-TL и Evolution-C - Инструкция по применению
ET	Evolution-TL and Evolution-C Removal Kits - Kasutusjuhend	UK	Набори екстракторів Evolution-TL і Evolution-C - Інструкція по застосуванню
FI	Evolution-TL- ja Evolution-C-irrotussarjat - Käyttöohjeet	KK	Evolution -TL мен Evolution-C тасымалды жанақар Қолданыс туралы нұсқау
EL	Kit αποίρεσης Evolution-TL και Evolution-C - Οδηγίες χρήσης	TK	Evolution-TL ve Evolution-C Çihaz Çıkartma Kitleri - Kullanım Talimatı
HU	Evolution-TL and Evolution-C Removal Kits - Instructions for Use	AR	Evolution-C و Evolution-TL أنظمة الإزالة تعليمات الاستخدام
LV	Evolution-TL un Evolution-C izņemšanas komplekti - Lietošanas instrukcija	TC	Evolution-TL Evolution-C 取出套件 使用說明書
LT	„Evolution-TL“ ir „Evolution-C“ stuburo implantų pašalinimo rinkiniai - Naudojimo instrukcijos	SC	Evolution-TL 和 Evolution-C 取出套件 使用说明书
NO	Evolution-TL un Evolution-C izņemšanas komplekti - Lietošanas instrukcija	KO	Evolution-TL 및 Evolution-C 제거 키트 사용 설명서
PL	Zestawy Evolution-TL Kit oraz Evolution-C Removal Kit-Instrukcja obsługi	ID	Evolution-TL dan Evolution-C - Petunjuk Penggunaan

Evolution-TL and Evolution-C Removal Kits

Instructions for Use

Manufacturer

Tecomet
5307 95th Avenue
Kenosha, WI 53144 USA

Description

The reusable manual orthopedic surgical instruments contained in the Evolution Removal Kits consist of bits for standard metric and English drives, hexalobe, broken and stripped screw removers, specific bits, extensions, ratchet handle, multi-axial handle, counter torque wrench, combination wrenches, needle-nose locking pliers and cutting instruments. The instruments are provided in an instrument storage container.

Intended Use

The Evolution-TL Removal Kit is intended for use in the removal of spinal hardware from the lumbar and thoracic regions of the spine. The Evolution-C Removal Kit is intended for use in the removal of spinal hardware from the cervical region of the spine.

Rx Only

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner.

Contraindications

These instruments are contraindicated in the removal of spinal hardware systems which use a pressure fit to retain its components or which twist off the threads of the locking screws.

Precautions

Instruments are provided non-sterile. Clean and sterilize prior to each use.

Disassemble instruments with removable parts prior to cleaning and sterilization.

Neutral pH cleaners are recommended.

Prior to use, care must be taken to protect nerves, vessels and/or organs from damage that may result from the use of these instruments.

Reprocessing Limitations

Repeated processing has minimal effect upon manual instruments. End of life is determined by wear and damage due to use.

Point of Use

Remove excess bodily fluids and tissue with a disposable wipe.

Instruments should not be allowed to dry prior to cleaning.

Containment / Transportation

Follow hospital protocols when handling contaminated and or biohazardous materials. Instrument should be cleaned within 30 minutes after use to minimize the potential of staining, damage, and drying after use.

Manual Cleaning Procedure

1. Prepare enzymatic neutral pH pre-soaks and cleaning agents according to manufacturer's instructions.
2. Rinse the devices with tap water for approximately two minutes while brushing with a soft bristled brush to remove most or all of the visible gross debris from the device. Careful attention should be paid to any pivots, threads, recesses, or crevices on the device. ***Note: Ratcheting T-handles must have their handles rotated 360 degrees forward and reverse after rinsing to insure that the instrument is fully drained.**
3. Place devices into an enzymatic detergent bath for five minutes. Scrub the devices using a soft bristled brush to remove any remaining debris. Careful attention should be paid to any pivots, threads, recesses, crevices or other hard to clean areas on the device. ***Note: Change the solution when it becomes grossly contaminated.**
4. Rinse devices for two additional minutes using tap water.
***Note: Ratcheting T-handles must have their handles rotated 360 degrees forward and reverse after rinsing to insure that the instrument is fully drained.**

5. Dry the instrument with a clean, disposable, absorbent, non-shedding wipe.

Automated Cleaning Procedure

A validation study for Automatic Cleaning has not been performed for these instruments. Please consult an appropriate responsible person or hospital protocol for direction.

Disinfection

Disinfection is only acceptable as an adjunct to full sterilization for reusable surgical instruments. See sterilization section below.

Inspection / Functional Testing

1. Carefully inspect each device to ensure that all visible blood and soil has been removed.
2. Visually inspect for damage and/or wear.
3. Check the action of moving parts (such as hinges and box-locks) to ensure smooth operation throughout the intended range of motion.
4. Check instruments with long slender features (particularly rotating instruments) for distortion.
5. Where instruments form part of a larger assembly, check that the devices assemble readily with mating components.
6. Inspect edges of the case, tray and cassette to ensure there aren't any sharp edges which might cut sterilization wrap. Inspect lid to assure there isn't any cracking, bending or warping. Verify lid continues to fit case tightly.

****Note: If damage or wear is noted that may compromise the function of the instrument, do not use the instrument and notify the appropriate responsible person.***

Sterilization

This procedure was validated to achieve a 10^{-6} Sterility Assurance Level (SAL). Because results may vary between autoclave equipment, it is highly recommended that the procedure be verified with the end-user's personal autoclave equipment using the standard test method described in the current ASTM standard (ASTM E 1766).

Instruments and the instrument case must be cleaned prior to sterilization. Parts should be arranged in a manner that they do not occlude the vent holes or the working end of devices and remain separate (not touching).

The case provided with the Evolution Removal Kits has been validated for loads under 19 lbs. **Do not overload the case.** The case, with instruments secured within the case, must be wrapped using CSR (sterilization wrap) following the AAMI double wrap method (ANSI/AAMI ST46), or equivalent. It is not recommended to stack the case during sterilization. Do not place case on side or upside down. The use of nonabsorbent tray liners can cause condensation to pool.

Method – Pre-Vacuum Steam Sterilization

Set parameters at ten minutes exposure at a temperature of 132°C (270°F) with a minimum of four (4) vacuum pulses prior to actual exposure.

Method – Gravity Displacement Steam Sterilization

Expose instruments for 30 minutes at a temperature of 132° C (270° F).

Handle sterilized kit following normal operating room procedures.

Allow case to cool thoroughly for a minimum of 15 minutes after sterilization cycle is complete.

Maintenance

Lubricate hinges, threads and other moving parts with a commercial water-based surgical grade instrument lubricant (such as instrument milk) to reduce friction and wear. Follow lubricant manufacturer's instructions.

Storage

Store sterile, packaged instruments in a designated, limited access area that is well ventilated and provides protection from dust, moisture, insects, vermin, and temperature/humidity extremes. Examine sterile instrument packages closely prior to opening to ensure that there has been no loss of package integrity.

Instructions for use


Cautions: The user of this instrument set acknowledges that implant manufacturers are constantly improving their designs, which could lead to a fit concern with custom drivers.

Do not use the driver(s) if the instrument does not properly fit into the implanted spinal hardware.

Prior to inserting the driver into the implant, verify that all foreign material is removed from the driving surface. Failure to do so could result in the damaging of implants and/or removal instruments.

Inspect instruments for damage or wear prior to use. Failure to do so could result in the damaging of implants or removal instruments. ***Note: If damage or wear is noted that may compromise the function of the instrument, do not use the instrument and notify the appropriate responsible person.**


BROKEN SCREW REMOVAL (Combined trephine and Bone Screw Remover)

1. Measure the outside diameter of broken screw and select trephine/ Broken Screw Remover that matches screw size.
2. Place trephine in handle. Turn counter-clockwise around broken screw. The trephine should clear around the screw as well as bite into screw.
3. Continue turning in counter-clockwise direction until screw is completely removed.

***Note: There may be bone in-growth on the screw preventing removal of the screw.**

Broken Screw Removal (Separate trephine and Bone Screw Remover)

1. Measure the outside diameter of the screw and select trephine that matches screw size.
2. Place trephine in handle of choice provided in the removal kit. Position trephine over broken screw. Note: The Trephine should clear around the screw. Approximately 5mm of screw should be exposed to properly seat the Trephine. To expose the broken screw, use a clockwise and counter-clockwise motion while applying downward pressure. Advance the trephine in this manner until the trephine reaches the broken screw. Note: During this process the teeth of the trephine may become clogged with debris. Remove trephine and clean debris from tool.
3. Select the bone screw remover that matches the diameter of the bone screw to be removed. Place bone screw remover into handle and position bone screw remover over the exposed screw. Use a counter-clockwise motion while applying downward pressure. The bone screw remover will bite into the screw. Continue turning in a counter-clockwise direction until the screw is completely removed.

***Note: There may be bone in-growth on the screw preventing removal of the screw.**

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