 Orthofix Inc.  
3451 Plano Parkway  
Lewisville, Texas 75056-9453 U.S.A.  
1-214-937-3199  
1-888-298-5700  
www.orthofix.com

Australian Sponsor  
Emergo Australia  
Level 20, Tower II  
Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

 Medical Device Safety Service  
(MDSS)  
Schiffgraben 41  
30175 Hannover  
Germany  
www.mdss.com

---

Device System Name:  
**InSWing™ Interspinous Spacer**

 Click directory below for desired language

<b>English</b>	<b>en</b>	<b>2-4</b>
<b>Español</b>	<b>es</b>	<b>5-7</b>
<b>Français</b>	<b>fr</b>	<b>8-10</b>
<b>Italiano</b>	<b>it</b>	<b>11-13</b>
<b>Português</b>	<b>pt</b>	<b>14-16</b>
<b>Deutsch</b>	<b>de</b>	<b>17-19</b>

Orthofix Inc.  
3451 Plano Parkway  
Lewisville, Texas 75056-9453 U.S.A.  
1-214-937-3199  
1-888-298-5700  
www.orthofix.com

Australian Sponsor  
Emergo Australia  
Level 20, Tower II  
Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

**EC REP** Medical Device Safety Service  
(MDSS)  
Schiffgraben 41  
30175 Hannover  
Germany  
www.mdss.com

## English en

### Device System Name:

## InSWing™ Interspinous Spacer

### Description:

The InSWing Interspinous Spacer is an implantable spinal construct that is manufactured from PEEK and is designed to fit between the spinous processes of the vertebrae. InSWing is implanted through a unilateral approach. The self-deploying and self-locking wings are designed to prevent the implant from pulling out. The implant limits extension and the flexion-limiting band (optional) is intended to improve the stability of the treated intervertebral lumbar segment.

The system is composed of the following components, each provided in an individual sterile package:

- InSWing Spacer implant 8 mm, Catalogue No. 17-0008
- InSWing Spacer implant 10 mm, Catalogue No. 17-0010
- InSWing Spacer implant 12 mm, Catalogue No. 17-0012
- InSWing Spacer implant 14 mm, Catalogue No. 17-0014
- InSWing Spacer implant 16 mm, Catalogue No. 17-0016

The following components are provided sterile and are packaged together as Catalogue No. 17-0001: InSWing band (qty. of 1), and clip (qty. of 2).

The InSWing spacers are made out of PEEK (polyetheretherketone) as described by ASTM F2026, with titanium alloy (Ti-6Al-4V) markers as described by ASTM F136. The InSWing band is made out of 100% PET (polyethylene terephthalate). The InSWing clip is made out of titanium alloy (Ti-6Al-4V) as described by ASTM F136.

The InSWing Interspinous Spacer is intended to be placed between the spinous processes of the vertebrae via a unilateral approach.

### Indications for Use:

The InSWing Interspinous Spacer is indicated to relieve leg and back pain due to symptomatic Lumbar Spinal Stenosis (LSS) from levels L1/L2 to L4/L5 with or without neurogenic intermittent claudication. The device can be used in combination with decompressive procedures providing that the spinous processes remained intact.

### Contraindications:

Contraindications include, but are not limited to:

1. Morbid obesity.
2. Mental illness.
3. Alcoholism or drug abuse.
4. Pregnancy.
5. Metal sensitivity/allergies.
6. Severe osteopenia/osteoporosis.
7. Patients unwilling or unable to follow post-operative care instructions.
8. Active local and systemic inflammation conditions.
9. Patient allergy to implant material or to its elements.
10. Patient's poor general condition that has a high risk of endangering health condition and/or loss of life as a result of surgical procedure.
11. Paget's disease.
12. Kyphosis requiring surgical correction.
13. Any circumstances not listed under the heading Indications for Use.
14. Other general or patient related contraindications against surgical implant treatment which may preclude potential benefit from implantation.

### Potential Adverse Events:

All of the possible adverse events associated with spinal surgery without instrumentation are possible. With instrumentation, a listing of possible adverse events includes, but is not limited to:

1. Neurological injury.
2. Vascular or visceral injury.
3. Fracture of the spinous process.
4. Early or late loosening of any or all of the components.
5. Incorrect positioning of the implant.
6. Disassembly, breakage and/or bending of any or all components.
7. Foreign body (allergic) reaction to implants, debris, corrosion products, including metallosis, straining, tumor formation, and/or auto-immune disease.
8. Pressure on the skin from component parts in patients with inadequate tissue coverage over the implant possibly causing skin penetration, irritation, and/or pain.
9. Post-operative change in spinal curvature, loss of correction, height, and/or reduction.
10. Infection.
11. Pain, discomfort, or abnormal sensations due to the presence of the device.
12. Chronic inflammation.

13. Hemorrhage.
14. Cessation of any potential growth of the operated portion of the spine.
15. Death.
16. Potential risks identified with the use of the device system may require additional surgery.

### Warnings and Precautions:

1. The safety and effectiveness of the interspinous spacer have been established only for spinal conditions with symptomatic LSS of the L1/L2-L4/L5 levels.
2. Potential risks identified with the use of this device system, which may require additional surgery, include: device component fracture, neurological injury, and vascular or visceral injury.
3. Single use only. Reuse of devices labeled as single-use (e.g. implants, drills, tacks, trial rods) could result in injury or reoperation due to breakage or infection.
4. Non-sterile: instruments are sold non-sterile, and therefore must be sterilized before use.
5. Do not use the InSWing implant past the expiration date noted on the package.
6. Do not use damaged products or opened packages.
7. Do not re-sterilize implants.
8. The implant should be placed as close as possible to the lamina without removing the supraspinous ligament. The interspinous ligament should be properly removed, as per the InSWing operative technique.
9. The implantation of the InSWing Interspinous Spacer System should be performed only by experienced spinal surgeons with specific training in the use of this implant because it is a technically demanding procedure presenting a risk of serious injury to the patient.
10. Do not use the components of this system with implants of other material composition or components from different manufacturers unless specifically stated.
11. All implants are intended for SINGLE USE ONLY. Any used implant should be discarded. Even though the device may appear undamaged, it may have small defects and internal stress patterns that may lead to fatigue failure.

### MRI Compatibility Information:

The InSWing Interspinous Spacer has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of InSWing Interspinous Spacer in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

### Cleaning:

InSWing Interspinous Spacer implants are provided STERILE. Do not use the implants if the package is opened or damaged or if the expiration date has passed. Please discard all opened and unused implants. Do not re-sterilize an opened and unused implant.

All instruments must be thoroughly cleaned after each use. Cleaning may be done using validated hospital methods or following the validated cleaning process described below.

None of the instruments in the system require disassembly prior to cleaning.

### From Point of Use:

Whenever possible, do not allow blood, debris or body fluids to dry on instruments. For best results and to prolong the life of the surgical instrument reprocess immediately after use.

1. Remove excess body fluids and tissue from instruments with a disposable, non-shedding wipe. Place instruments in a basin of purified water or in a tray covered with damp towels. Do not allow saline, blood, body fluids, tissue, bone fragments or other organic debris to dry on instruments prior to cleaning.
2. For optimal results, instruments should be cleaned within 30 minutes of use or after removal from solution to minimize the potential for drying prior to cleaning.
3. Used instruments must be transported to the central supply in closed or covered containers to prevent unnecessary contamination risk.

**Note:** Soaking in proteolytic enzymatic detergents or other pre-cleaning solutions facilitates cleaning, especially in instruments with complex features and hard-to-reach areas (e.g. cannulated and tubular designs, etc.). These enzymatic detergents as well as enzymatic foam sprays break down protein matter and prevent blood and protein based materials from drying on instruments. Manufacturer's instructions for preparation and use of these solutions should be explicitly followed.

### Preparation for Cleaning:

1. All instruments with moving parts (e.g., knobs, triggers, hinges) should be placed in the open position to allow access of the cleaning fluid to areas that are difficult to clean.
2. Soak the instruments for a minimum of 10 minutes in purified water prior to the manual or automated cleaning process.
3. Use a soft cloth or a soft plastic bristle brush to remove any visible soil from the instruments prior to manual or automated cleaning. Use a soft plastic bristle brush or a pipe cleaner to remove soil from any inner lumens. You can also use a syringe (if appropriate) for hard to reach areas.

- Enzymatic detergent should be used for manual and automated cleaning. All enzymatic detergents should be prepared at the use dilution and temperature recommended by the manufacturer. Softened tap water may be used to prepare the enzymatic detergents. Use of recommended temperatures is important for optimal performance of enzymatic detergent.

#### Manual Cleaning:

- Completely submerge instruments in an enzymatic detergent and allow to soak for 20 minutes. Use a soft-bristled, nylon brush to gently scrub the device until all visible soil has been removed. Particular attention must be given to crevices, lumens, mated surfaces, connectors and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft-bristled brush (i.e. pipe cleaner brush).
- Remove the instruments from the enzymatic detergent and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes and other difficult to reach areas.
- Place prepared cleaning solution in a sonication unit. Completely submerge device in cleaning solution and sonicate for 10 minutes.
- Rinse instrument in purified water for at least 3 minutes or until there is no sign of blood or soil on the device or in the rinse stream. Thoroughly and aggressively flush lumens, holes and other difficult to reach areas.
- Repeat the sonication and rinse steps above.
- Remove excess moisture from the instrument with a clean, absorbent and non-shedding wipe.
- Inspect the instruments for visible soil.
- If visible soil is noted, repeat the steps listed above.

#### Automated Cleaning:

- Completely submerge the instruments in an enzymatic detergent and allow to soak and sonicate for 10 minutes each. Use a soft nylon bristled brush to gently scrub the device until all visible soil has been removed. Particular attention must be given to crevices, lumens, mated surfaces, connectors and other hard to clean areas. Lumens should be cleaned with a long, narrow, soft nylon bristled brush (i.e. pipe cleaner). Use of a syringe or water jet will improve flushing of difficult to reach areas and closely mated surfaces.
- Remove instruments from the cleaning solution and rinse in purified water for a minimum of 1 minute. Thoroughly and aggressively flush lumens, blind holes and other difficult to reach areas.
- Place instruments in a suitable washer/disinfectant basket and process through a standard instrument washer/disinfectant cleaning cycle.
- Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
- The following minimum parameters are essential for thorough cleaning.
  - 2 minute prewash with cold tap water
  - 1 minute prewash with hot tap water
  - 2 minute detergent wash with hot tap water (64-66°C/146-150°F)
  - 1 minute hot tap water rinse
  - 2 minute thermal rinse with purified water (80-93°C/176-200°F)
  - 1 minute purified water rinse (64-66°C/146-150°F)
  - 7 to 30 minute hot air dry (116°C/240°F)
- Inspect the instruments for visible soil.
- If visible soil is noted, repeat the above listed steps until no visible soil is noted.

**Note:** Certain cleaning solutions such as those containing caustic soda, formalin, glutaraldehyde, bleach, and/or other alkaline cleaners may damage instruments. These solutions should not be used.

**Note:** Visually inspect instruments after cleaning and prior to each use. Discard or return to Orthofix any instruments that are broken, discolored, corroded, have cracked components, pits, gouges, or are otherwise found defective. Do not use defective instruments.

#### Torque Limiting Instrument Maintenance

- If a torque-limiting handle has been dropped, impacted or used incorrectly, return to Orthofix.
- Torque-limiting handles require maintenance at minimum, every three years or per your service agreement. Please return your torque limiting handles to Orthofix for required maintenance.

#### Instrument End of Life Determination:

Do not reuse Single Use instruments. Visually inspect the reusable instruments to determine if the instrument has reached end of life. Orthofix reusable instruments have reached End of Life when:

- Instruments show signs of damage such as binding, bending, breakage, overt signs of wear and/or any other conditions which may impact the devices safe and effective use.
- Instruments intended for cutting bone and/or tissue (e.g. tap, rasp, curette, rongeur) – when any of the cutting surfaces show signs of wear such as nicks, abrasions or otherwise dulled cutting surfaces.
- Instruments that interface with other devices (e.g. implants, instruments, handles) - when the mating feature binds, fails to attach or fails to hold the device securely. The instrument function should be verified prior to each use.
- Do not use instruments which reached End of Life. Discard End of Life instruments per your hospital procedure or return to Orthofix for disposal.

#### Sterilization:

The InSWing Interspinous Spacer instruments (medical-grade stainless steel and silicone handles) should be steam sterilized by the hospital using the following recommended cycles when utilizing a sterilization wrap:

Method: Steam	or:	Method: Steam
Cycle: Gravity		Cycle: Prevac
Temperature: 250°F (121°C)		Temperature: 270°F (132°C)
Exposure time: 30 minutes		Exposure time: 8 minutes

The InSWing Interspinous Spacer implants are sterilized using gamma irradiation sterilization. Do not re-sterilize.

#### Physician Information:

##### Preoperative:

- Only patients that meet the criteria described in the indications should be selected.
- Patient conditions and/or predispositions such as those addressed in the aforementioned contraindications should be avoided.
- Correct selection of the implant is extremely important.
- Use care in handling and storage of implant components. Inspection should be made to determine if components have been damaged during storage or previous procedures.
- An adequate inventory of implant sizes should be available at the time of surgery.
- Certain special surgical instruments are required to perform this surgery. Review of the use and handling of these instruments is very important.

##### Intraoperative:

- A unilateral surgical approach is the recommended access to the interspinous space in order to minimize trauma to the soft tissue.
- Use great care to insure that the spinous process scraper is used to clean the remains of the interspinous ligament from spinous process.
- Do not damage the supraspinous ligament. Be sure to detach muscle on the contralateral side with the scraper to allow the wings to properly deploy
- Correct sizing is imperative for good clinical outcome. When using the sizing trial or distractor, be careful not to over distract. Too much distraction can induce kyphosis.
- During implantation of the InSWing, do not use excessive force as it may cause a fracture of the spinous process.
- Before closing, use fluoroscopy to verify the final implant position. Complete the surgery using standard posterior lumbar closure procedures.

##### Postoperative:

- Detailed instructions on the use and limitations of the implant should be given to the patient. The patient must be made aware of the limitations of the implant.
- Periodic X-rays for at least the first year postoperatively are recommended for close comparison with postoperative conditions to detect any evidence of changes in position, bending or cracking of components.
- Surgical implants must never be reused. Any retrieved devices should never be reused in another surgical procedure. The retrieved parts should be handled and disposed of in such a manner as to ensure that reuse is not possible.

##### Removal of Implant:

In certain circumstances, removal may be necessary. These reasons include, but are not limited to:

- Dismantling or breakage of implant.
- Fracture of spinous process with pain.
- Infection.
- Surgery involving spinous processes.

The recommended surgical steps for removal are as follows:

- Cut the band and knot (if present).
- Remove the band and clip used to secure the original knot.
- Make a small opening in the intact fascia on the side of the opening wings.
- Insert a raspator (periosteal elevator) along the spinous process and close (lift) one of the wings of the implant while pulling the implant out. This can also be done with a finger.

##### Packaging:

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The InSWing Interspinous Spacer implants and instruments should be stored in a cool and dry place.

The InSWing Interspinous Spacer instruments are provided in a modular case specifically intended to contain and organize the system's components. The system's instruments are organized into trays within each modular case for easy retrieval during surgery. These trays also provide protection to the system components during shipping. Additionally, individual instruments are provided in sealed poly bags with individual product labels.

The InSWing Interspinous Spacer implants are supplied sterile and packaged in a dual barrier system to facilitate introduction into the sterile field in the surgical theater. The devices are considered sterile, unless the sterile expiration date is elapsed or the package is opened or damaged, e.g., loss of seal integrity or is not in the original packaging.

##### Product Complaints:

Any Health Care Professional (e.g., customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Orthofix Inc., 3451 Plano Parkway, Lewisville, TX 75056, USA, by telephone at 1-214-937-3199 or 1-888-298-5700 or by e-mail at complaints@orthofix.com.


##### Further Information:

A recommended operative technique for the use of this system is available upon request from Orthofix at the phone numbers provided above.

##### Latex Information:

The implants, instruments and/or packaging material for the InSWing Interspinous Spacer System are not formulated with and do not contain natural rubber. The term "natural rubber" includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

**Caution:** The InSWing Interspinous Spacer has not been cleared by the US Food and Drug Administration and is not available in the United States.

	See Instructions for Use		Use By Date
	Orthofix.com/IFU		Manufacturer
	Single Use Only Do Not Reuse		Authorized Representative
	Catalogue Number		Do Not Resterilize
	Lot Number		Serial Number
	Sterilized Using Irradiation		Provided Non-Sterile

Orthofix Inc.  
3451 Plano Parkway  
Lewisville, Texas 75056-9453 U.S.A.  
1-214-937-3199  
1-888-298-5700  
www.orthofix.com

Australian Sponsor  
Emergo Australia  
Level 20, Tower II  
Darling Park  
201 Sussex Street  
Sydney, NSW 2000  
Australia

**EC REP** Medical Device Safety Service  
(MDSS)  
Schiffgraben 41  
30175 Hannover  
Germany  
www.mdss.com

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