Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

Australian Sponsor
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

 Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

Device System Name:
SambaScrew® SI Fixation System

Click directory below for desired language

English  **2-3**

Español  **4-6**

Français  **7-9**

Italiano  **10-12**

Deutsch  **13-15**

Português  **16-18**

Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

Australian Sponsor
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

EC REP Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

English EN

Device System Name:

SambaScrew® SI Fixation System

Description:

The SambaScrew SI Fixation System is a temporary, multiple component system consisting of non-sterile instruments and non-sterile, medical-grade titanium (Ti-6Al-4V ELI) implants. The system consists of a 9mm diameter, cannulated screw with multiple orifices on its shaft. The Steinmann Pins, Variable Drill Bit, and Packing Tube are single use devices and should be discarded after use.

Indications for Use:

The SambaScrew SI Fixation System is intended for fixation of sacroiliac joint disruptions. This fixation device is to only be used in skeletally mature patients.

Contraindications:

The SambaScrew SI Fixation System is contraindicated for use in patients with:

1. Open wounds, infection, presence of tumor, pregnancy, osteoporosis, certain metabolic disorders affecting osteogenesis, certain inflammatory/neuromuscular conditions, and certain neuromuscular deficits which would place an unusually heavy load on the device during the healing period.
2. The implant is made from Ti-6Al-4V ELI (medical-grade titanium alloy). The fixation implant is contraindicated in any individual with a known or suspected allergy, sensitivity or intolerance to metal.

Potential Adverse Events:

Potential adverse events include, but are not limited to:

1. Allergic reaction or metal sensitivity to foreign body.
2. Cardiovascular system compromise.
3. Death.
4. Decrease in bone density due to stress shielding.
5. Device bending, disassembly, fracture, loosening, migration and/or retropulsion, or subsidence.
6. Dural tears, neural structure injury.
7. Fixation implant migration with or without bone fracture.
8. Fracture of pelvis or sacrum.
9. Gastrointestinal complications (i.e., ileus or bowel perforation).
10. Hemorrhage.
11. Incisional complications (i.e., dehiscence, hematoma).
12. Infection (Incisional or implant site).
13. Loss of spinal mobility or function.
14. Loosening or fracture of fixation implant.
15. Malfunction of fixation device and/or instruments.
16. Malposition of the fixation device.
17. Migration of fixation implant.
18. Neurological injury/deficit which may range from paresthesias to muscle paralysis, loss of rectal or bladder sphincter control, radiculopathies.
19. Organ, connective tissue or nerve damage.
20. Osteoarthritis.
21. Pain, discomfort or abnormal sensation due to device presence.
22. Persistent low back pain.
23. Reproductive system compromise.
24. Screw back-out or breakage possibly leading to local pain, perforation or irritation of adjacent structures.
25. Sepsis.
26. Urological compromise (i.e., infection/retention).
27. Vascular injury.
28. Wound hematoma.

Note: Additional surgery may be necessary to correct some of these anticipated adverse events.

Warnings and Precautions:

1. The SambaScrew device should only be used by surgeons who have been trained in the use of this device. Information on laboratory and clinical training, as well as additional brochures with a detailed description of proper surgical technique, may be obtained from Orthofix. See the SambaScrew SI Fixation System Surgical Manual for instructions on the implant procedure.
2. Infection may occur immediately following implant fixation or a long time afterwards due to transient bacteremia such as caused by dental treatment(s), endoscopic examination or any other minor surgical procedure. To avoid infection at the implant fixation site, it may be advisable to use antibiotic prophylaxis before and/or after such procedures.
3. Single use only. Reuse of devices labeled as single-use (e.g. implants, drills, tacks, trial rods) could result in injury or reoperation due to breakage or infection.
4. All implants are intended for SINGLE USE ONLY. Any used implant should be discarded. Even though the device may appear undamaged, it may have small defects and internal stress patterns that may lead to fatigue failure.

MRI Compatibility Information:

The SambaScrew SI Fixation System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of SambaScrew SI Fixation System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Cleaning:

SambaScrew SI Fixation System instruments and implants are provided clean but not sterile. Once an implant comes in contact with any human tissue or bodily fluid, it should not be re-sterilized or used. Please discard all contaminated implants.

All instruments must be thoroughly cleaned after each use. Cleaning may be done using validated hospital methods or following the validated cleaning processes described below.

From Point of Use:

Whenever possible, do not allow blood, debris or body fluids to dry on instruments. For best results and to prolong the life of the surgical instrument reprocess immediately after use.

1. Remove excess body fluids and tissue from instruments with a disposable, non-shedding wipe. Place instruments in a basin of purified water or in a tray covered with damp towels. Do not allow saline, blood, body fluids, tissue, bone fragments or other organic debris to dry on instruments prior to cleaning.
2. For optimal results, instruments should be cleaned within 30 minutes of use or after removal from solution to minimize the potential for drying prior to cleaning.
3. Used instruments must be transported to the central supply in closed or covered containers to prevent unnecessary contamination risk.

Note: Soaking in proteolytic enzymatic detergents or other pre-cleaning solutions facilitates cleaning, especially in instruments with complex features and hard-to-reach areas (e.g. cannulated and tubular designs, etc.). These enzymatic detergents as well as enzymatic foam sprays break down protein matter and prevent blood and protein based materials from drying on instruments. Manufacturer's instructions for preparation and use of these solutions should be explicitly followed.

Preparation for Cleaning:

1. All instruments with moving parts (e.g., knobs, triggers, hinges) should be placed in the open position to allow access of the cleaning fluid to areas that are difficult to clean.
2. Soak the instruments for a minimum of 10 minutes in purified water prior to the manual or automated cleaning process.
3. Use a soft cloth or a soft plastic bristle brush to remove any visible soil from the instruments prior to manual or automated cleaning. Use a soft plastic bristle brush or a pipe cleaner to remove soil from any inner lumens. You can also use a syringe (if appropriate) for hard to reach areas.
4. Enzymatic detergent should be used for manual and automated cleaning. All enzymatic detergents should be prepared at the use dilution and temperature recommended by the manufacturer. Softened tap water may be used to prepare the enzymatic detergents. Use of recommended temperatures is important for optimal performance of enzymatic detergent.

The following SambaScrew SI Fixation System instruments requires disassembly prior to cleaning

Instructions for Disassembly:

1. Parallel Guide Set (PN 11-1109-FD7):
There are three parts to the Parallel Guide Set: 1. Locking Nut, 2. Sliding Block, and 3. Parallel Guide Frame. The Parallel Guide Set requires disassembly prior to cleaning using the following steps:
a. Unthread the Nut from the Sliding Block.
b. Remove the Sliding Block from the Parallel Guide Frame.
2. Packing Plunger (11-1109-FD8-1) and Packing Tube (11-1109-FD8-2):
There are two parts: 1. Packing Plunger and 2. Packing Tube (single use only). The Packing Plunger and Packing Tube require disassembly prior to cleaning using the following steps:
a. Remove the Plunger from the Tube.
b. Discard the Tube, which is single use only. Do not discard the Plunger.
3. Variable Drill Bit (PN 11-1109-FD9-1) and Adjustable Drill Collar (PN 11-1109-FD9-2):
There are two parts: 1. Variable Drill Bit (single use only) and 2. Drill Collar. The Variable Drill Bit and Adjustable Drill Collar require disassembly prior to cleaning using the following steps:
a. Disengage the Drill Collar from the Variable Drill Bit's shaft.
b. Slide the Drill Collar over the Drill Bit Shank to take it off.
c. Discard the Variable Drill Bit, which is single use only. Do not discard the Adjustable Drill Collar.

Manual Cleaning:

1. Completely submerge instruments in an enzymatic detergent and allow to soak for 20 minutes. Use a soft-bristled, nylon brush to gently scrub the device until all visible soil has been removed. Particular attention must be given to crevices, lumens, mated surfaces, connectors and other hard-to-clean areas. Lumens should be cleaned with a long, narrow, soft-bristled brush (i.e. pipe cleaner brush).
2. Remove the instruments from the enzymatic detergent and rinse in tap water for a minimum of 3 minutes. Thoroughly and aggressively flush lumens, holes and other difficult to reach areas.
3. Place prepared cleaning solution in a sonication unit. Completely submerge device in cleaning solution and sonicate for 10 minutes.
4. Rinse instrument in purified water for at least 3 minutes or until there is no sign of blood or soil on the device or in the rinse stream. Thoroughly and aggressively flush lumens, holes and other difficult to reach areas.
5. Repeat the sonication and rinse steps above.
6. Remove excess moisture from the instrument with a clean, absorbent and non-shedding wipe.
7. Inspect the instruments for visible soil.
8. If visible soil is noted, repeat the steps listed above.

Automated Cleaning:

1. Completely submerge the instruments in an enzymatic detergent and allow to soak and sonicate for 10 minutes each. Use a soft nylon bristled brush to gently scrub the device until all visible soil has been removed. Particular attention must be given to crevices, lumens, mated surfaces, connectors and other hard to clean areas. Lumens should be cleaned with a long, narrow, soft nylon bristled brush (i.e. pipe cleaner). Use of a syringe or water jet will improve flushing of difficult to reach areas and closely mated surfaces.
2. Remove instruments from the cleaning solution and rinse in purified water for a minimum of 1 minute. Thoroughly and aggressively flush lumens, blind holes and other difficult to reach areas.
3. Place instruments in a suitable washer/disinfector basket and process through a standard instrument washer/disinfector cleaning cycle.
4. Orient instruments into the automated washer's carriers as recommended by the washer manufacturer.
5. The following minimum parameters are essential for thorough cleaning.
 - a. 2 minute prewash with cold tap water
 - b. 1 minute prewash with hot tap water
 - c. 2 minute detergent wash with hot tap water (64-66°C/146-150°F)
 - d. 1 minute hot tap water rinse
 - e. 2 minute thermal rinse with purified water (80-93°C/176-200°F)
 - f. 1 minute purified water rinse (64-66°C/146-150°F)
 - g. 7 to 30 minute hot air dry (116°C/240°F)
6. Inspect the instruments for visible soil.
7. If visible soil is noted, repeat the above listed steps until no visible soil is noted.

Note: Certain cleaning solutions such as those containing caustic soda, formalin, glutaraldehyde, bleach, and/or other alkaline cleaners may damage instruments. These solutions should not be used.

Instructions for Assembly and Inspection:

1. Parallel Guide Set (PN 11-1109-FD7):
There are three parts to the Parallel Guide Set: 1. Locking Nut, 2. Sliding Block, and 3. Parallel Guide Frame. The Parallel Guide Set requires assembly after cleaning using the following steps:
 - a. Slide the Sliding Block through the underside of the Parallel Guide Frame with the threads leading.
 - b. Slide the Nut over the Sliding Block's shaft and thread it into place, making sure the Sliding Block seats properly.
2. Packing Plunger (PN 11-1109-FD8-1) and Packing Tube (PN 11-1109-FD8-2):
There are two parts: 1. Packing Plunger and 2. Packing Tube (single use only). The Packing Plunger and Packing Tube require assembly after cleaning using the following step:
 - a. Insert the Plunger into the Tube.
3. Variable Drill Bit (PN 11-1109-FD9-1) and Adjustable Drill Collar (PN 11-1109-FD9-2):
There are two parts: 1. Variable Drill Bit (single use only) and 2. Adjustable Drill Collar. The Variable Drill Bit and Adjustable Drill Collar require assembly after cleaning using the following steps:
 - a. Slide the Drill Collar, with the arrow pointing towards the proximal end of the Drill Bit, over the Drill Bit Shank until the desired depth position is engaged.
 - b. Flex the parallel beams to release the Drill Collar from the detent engagement.

Note: Visually inspect instruments after cleaning and prior to each use. Discard or return to Orthofix any instruments that are broken, discolored, corroded, have cracked components, pits, gouges, or are otherwise found defective. Do not use defective instruments.

Instrument End of Life Determination:

Do not reuse Single Use instruments. Visually inspect the reusable instruments to determine if the instrument has reached end of life. Orthofix reusable instruments have reached End of Life when:

1. Instruments show signs of damage such as binding, bending, breakage, overt signs of wear and/or any other conditions which may impact the devices safe and effective use.
2. Instruments intended for cutting bone and/or tissue (e.g. tap, rasp, curette, rongeur) – when any of the cutting surfaces show signs of wear such as nicks, abrasions or otherwise dulled cutting surfaces.
3. Instruments that interface with other devices (e.g. implants, instruments, handles) - when the mating feature binds, fails to attach or fails to hold the device securely. The instrument function should be verified prior to each use.
4. Do not use instruments which reached End of Life. Discard End of Life instruments per your hospital procedure or return to Orthofix for disposal.

Sterilization:

The SambaScrew SI Fixation System implants and instruments are provided NON-STERILE and must be thoroughly cleaned and sterilized prior to each use.

AW-11-9901 Rev.AC
SA-1904-PL-US © Orthofix, Inc. 4/2020

Place the dry instruments in the designated locations in the case (see markings on tray surface that indicate proper storage location in the case).

Sterilization for Instrument Cases with Blue Wrap:

It is recommended to double wrap instrument cases in standard central supply FDA cleared wrap prior to steam sterilization using the gravity cycles. A single wrap (or pouch) may be used for the pre-vacuum cycle.

Only sterile products should be placed in the operative field. For a 10⁻⁶ sterility assurance level, these products are recommended to be steam sterilized by the hospital using one of the three sets of process parameters below:

Sterilization Parameters for Use in the United States for Instrument Cases:

Method: Steam
Cycle: Pre-Vacuum
Temperature: 270°F (132°C)
Exposure time: 4 Minutes
Drying time: NA – Immediate Use Cycle

The following gravity sterilization cycles are not considered by the United States Food and Drug Administration (US FDA) to be standard sterilization cycles. Users should only use sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization containers) that have been cleared by the US FDA for the following selected gravity sterilization cycle specification.

Method: Steam*	or:	Method: Steam
Cycle: Gravity		Cycle: Gravity
Temperature: 273°F (134°C)		Temperature: 250°F (121°C)
Exposure time: 20 Minutes*		Exposure time: 60 Minutes
Drying time: 15-30 Minutes		Drying time: 15-30 Minutes

Note: Because of the many variables involved in sterilization, each medical facility should calibrate and verify the sterilization process (e.g., temperatures and times) used for their equipment.

*Some non-U.S. Health Care Authorities recommend sterilization according to these parameters in order to minimize the potential risk of transmission of Creutzfeldt-Jakob disease, especially for surgical instruments that could come into contact with the central nervous system. This sterilization cycle has not been validated for prion inactivation.

Sterilization for Screw Caddy with Blue Wrap:

The SambaScrew SI Fixation System implants are supplied NON-STERILE. Prior to use, all implants should be placed in the appropriate Orthofix caddy which will be wrapped in a FDA cleared sterilization wrap and placed in the autoclave for sterilization by the hospital using one of the following recommended cycles:

Method: Steam	or:	Method: Steam
Cycle: Gravity		Cycle: Pre-vacuum (minimum 4 pulses)
Temperature: 270°F (132°C)		Temperature: 270°F (132°C)
Exposure time: 15 minutes		Exposure time: 4 minutes
Drying time: 30 minutes		Drying time: 30 minutes

Validation and routine monitoring should be performed per ANSI/AAMI ST79 Comprehensive guide to steam sterilization and sterility assurance in health care facilities. Other cycles may be used as long as they comply with the above practices and provide a sterility assurance level of 10⁻⁶.

Packaging:

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The SambaScrew SI Fixation System instruments and implants are provided in a modular case and caddy specifically intended to contain and organize the system's components. The system's instruments are organized in a modular case for easy retrieval during surgery. The modular case also provides protection to the system components during shipping. Additionally, individual instruments and implants are provided in sealed poly bags with individual product labels.

Product Complaints:

Any Healthcare Professional (e.g., customer or user of this system) who has any complaints or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Orthofix Inc., 3451 Plano Parkway, Lewisville, TX 75056, USA, Telephone: 214-937-3199 or 1-888-298-5700 or via email at complaints@orthofix.com.

Further Information:

A recommended operative technique for the use of this system is available upon request from Orthofix at the phone numbers provided above.

Latex Information:

The implants, instruments and/or packaging material for the SambaScrew SI Fixation System are not formulated with and do not contain natural rubber. The term "natural rubber" includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

Caution: Federal law (USA) restricts these devices to sale by or on the order of a physician.

Rx Only		Federal (U.S.A.) law restricts this device to sale by or on the order of a physician	
	See Instructions for Use		Manufacturer
	Single Use Only Do Not Reuse		Authorised Representative
	Catalogue Number		Lot Number
	Provided Non-Sterile		

 Orthofix Inc.
3451 Plano Parkway
Lewisville, Texas 75056-9453 U.S.A.
1-214-937-3199
1-888-298-5700
www.orthofix.com

Australian Sponsor
Emergo Australia
Level 20, Tower II
Darling Park
201 Sussex Street
Sydney, NSW 2000
Australia

 Medical Device Safety Service
(MDSS)
Schiffgraben 41
30175 Hannover
Germany
+49 511 6262 8630
www.mdss.com

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