**INSTRUCTIONS FOR USE**  
Important Information – Please Read Prior to Use

**Rx Only**  
**CE** 0086

 Orthofix Inc.  
3451 Plano Parkway  
Lewisville, Texas 75056-9453 U.S.A.  
1-214-937-3199  
1-888-298-5700  
www.orthofix.com

**EC REP** Medical Device Safety Service  
(MDSS)  
Schiffgraben 41  
30175 Hannover  
Germany  
+49 511 6262 8630  
www.mdss.com

---

Device System Name:  
**NGage<sup>®</sup> Surgical Mesh System**

 Click directory below for desired language

<b>English</b> <b>EN</b>	<b>2-3</b>
<b>Español</b> <b>ES</b>	<b>4-5</b>
<b>Italiano</b> <b>IT</b>	<b>6-7</b>
<b>Deutsch</b> <b>DE</b>	<b>8-9</b>
<b>Português</b> <b>PT</b>	<b>10-11</b>
<b>英文</b> <b>ZH-TW</b>	<b>12-13</b>

# INSTRUCTIONS FOR USE

Important Information – Please Read Prior to Use

Rx Only  
CE 0086

 Orthofix Inc.  
3451 Plano Parkway  
Lewisville, Texas 75056-9453 U.S.A.  
1-214-937-3199  
1-888-298-5700  
www.orthofix.com

 Medical Device Safety Service  
(MDSS)  
Schiffgraben 41  
30175 Hannover  
Germany  
+49 511 6262 8630  
www.mdss.com

## English EN

### Device System Name:

## NGage® Surgical Mesh System

#### Description:

The NGage Surgical Mesh System is a diamond pattern, surgical mesh device. The body is manufactured from commercially pure (CP) titanium conforming to ASTM F67, the end rings, the standard ring, and the screws are manufactured from titanium alloy that conforms to ASTM F136. Because of the construction, the angle and the length of the mesh can be reduced incrementally to adjust it to individual anatomical conditions. The system is sold non-sterile.

#### Indications for Use:

The NGage Surgical Mesh System is indicated for use in the thoraco-lumbar spine (T1-L5) to replace a diseased vertebral body resected or excised for the treatment of tumors, to achieve anterior decompression of the spinal cord and neural tissues, and to restore the height of a collapsed vertebral body.

The NGage Surgical Mesh System is also indicated for treating fractures of the thoracic and lumbar spine.

The NGage Surgical Mesh System is designed to restore the biomechanical integrity of the anterior, middle, and posterior spinal column even in the absence of fusion for a prolonged period. It is recommended to pack bone graft material inside the mesh cage prior to implantation.

The NGage Surgical Mesh System is intended for use with supplemental fixation. As an example, the supplemental fixation system that may be used with the system is the Orthofix Spinal Fixation System.

#### Contraindications:

The NGage Surgical Mesh System, as with other metallic orthopedic appliances, is contraindicated for use in patients with active infections in which the use of an implant could preclude adequate and appropriate treatment of the infection. The system is also contraindicated for use in patients with known or suspected metal allergies.

#### Potential Adverse Events:

Potential adverse events include, but are not limited to:

1. Failure of the device to provide adequate mechanical stability.
2. Loss of fixation of the implant.
3. Device component failure.
4. Migration or bending of the device.
5. Loss of bony alignment.
6. Non-union.
7. Fracture of bony structures.
8. Resorption without incorporation of any bone graft utilized.
9. Immunogenic response to the implant materials.

**Note:** As with any surgical procedure, there are risks involved in orthopedic surgery. Infrequent operative and postoperative complications known to occur are: early or late infection which may result in the need for additional surgeries, damage to blood vessels, spinal cord or peripheral nerves, pulmonary emboli, loss of sensory and/or motor function, impotence, and permanent pain and/or deformity. Rarely, some complications may be fatal.

#### Warnings and Precautions:

The surgeon should be aware of the following when using metallic implants:

1. The correct selection of the implant is extremely important. The potential for success is increased by the selection of the proper size, shape and design of the implant. The size and shape of human bones present limiting restrictions on the size and strength of implants. No implant can be expected to withstand the unsupported stresses of full weight bearing.
2. Correct handling of the implant is extremely important. Contouring the metal devices is to be avoided where possible. If contouring is necessary, it should not be bent sharply, reverse bent, notched or scratched. All of these operations can produce defects in surface finish and internal stress concentrations, which may become the focal point for eventual failure of the device.
3. Single use only. No metallic surgical implant should be reused. Any metal implant once used should be discarded. Even though the device appears undamaged, it may already have small defects and internal stress patterns which may lead to fatigue failure.
4. Non-sterile; the surgical mesh and instruments are sold non-sterile, and therefore, must be sterilized before each use.
5. Postoperative care is important. The patient should be instructed in the limitations of the metallic implant and should be cautioned regarding weight bearing and body stress on the device prior to secure bone healing.
6. Reuse of devices labeled as single-use could result in injury or re-operation due to breakage or infection. Do not attempt to re-sterilize single-use implants that come in contact with body fluids.

#### MRI Compatibility Information:

The NGage Surgical Mesh System has not been evaluated for safety and compatibility in the Magnetic Resonance (MR) environment. This system has not been tested for heating or migration in the MR environment.

#### Cleaning:

All instruments and implants must first be cleaned using established hospital methods before sterilization and introduction into a sterile field. Additionally, all instruments and implants that have been previously taken into a sterile surgical field must first be cleaned using established hospital methods before sterilization and reintroduction into a sterile surgical field. Cleaning can include the use of neutral cleaners followed by a deionized water rinse. All products should be treated with care. Improper use or handling may lead to damage and possible improper functioning of the device.

#### Sterilization:

The NGage Surgical Mesh System components are supplied NON-STERILE. Prior to use, all components should be steam sterilized by the hospital using the recommended cycle:

Method: Steam	or:	Method: Steam
Cycle: Gravity		Cycle: Prevac
Temperature: 250°F (121°C)		Temperature: 270°F (132°C)
Exposure time: 30 minutes		Exposure time: 8 minutes

**Packaging:**

Packages for each of the components should be intact upon receipt. If a consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked for damage prior to use. Damaged packages or products should not be used and should be returned to Orthofix.

The NGage Surgical Mesh System instruments and implants are provided in modular cases specifically intended to contain and organize the system's components. The system's instruments are organized into trays within each modular case for easy retrieval during surgery. These trays also provide protection to the system components during shipping. Additionally, individual instruments and implants are provided in sealed poly bags with individual product labels.

**Product Complaints:**

Any Healthcare Professional (e.g., customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction with the product quality, identity, durability, reliability, safety, effectiveness and/or performance, should notify Orthofix Inc., 3451 Plano Parkway, Lewisville, TX 75056, USA, by telephone at 1-214-937-3199 or 1-888-298-5700 or by e-mail at [complaints@orthofix.com](mailto:complaints@orthofix.com).

**Further Information:**

A recommended operative technique for the use of this system is available upon request from Orthofix at the phone numbers provided above.

**Latex Information:**

The implants, instruments and/or packaging material for the NGage Surgical Mesh System are not formulated with and do not contain natural rubber. The term "natural rubber" includes natural rubber latex, dry natural rubber, and synthetic latex or synthetic rubber that contains natural rubber in its formulation.

**Caution:** Federal law (USA) restricts these devices to sale by or on the order of a physician.

<b>Rx Only</b>	Federal (U.S.A.) law restricts this device to sale by or on the order of a physician	
	See Instructions for Use	 Manufacturer
	Orthofix.com/IFU	 Authorized Representative
	Single Use Only Do Not Reuse	 Lot Number
	Provided Non-Sterile	 Catalogue Number

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