

FORZA[™]
SPACER SYSTEM




FORZA™

S P A C E R S Y S T E M

The FORZA™ Spacer System has been designed to help optimize Transforaminal Lumbar Interbody Fusion (TLIF) and Posterior Lumbar Interbody Fusion (PLIF) procedures with surgeon-designed implants and instruments.


Bullethead Spacers:

- Bulleted nose for assisting with distraction.
- Threaded connection to the Implant Inserter with grooves for a secure hold that is firm enough for an insert and rotate technique.
- Chamfered edges for easier rotation with the straight cage.
- Narrow double bullet nose on the straight cage regardless of the size to assist with distraction.
- Aggressive anti-migration ribs for resisting implant migration.
- Full height vertical tantalum markers 1mm from the ends for clear in-situ implant positioning.
- Provided sterile for inventory efficiency.




9mm Straight Spacer 0°

11mm Straight Spacer 8°




9mm Curved Spacer 0°

11mm Curved Spacer 8°



Streamlined Instrumentation:


- Firm connection to the implant, even for an insert and rotate technique with the straight cage.
- All Inserters provide great visualization capability.
- Modular Implant Inserter with different handle angle locations to match surgeon preference.
- MIS Inserter for improved visualization in tight spaces.
- Monolithic implant trials for speed and efficiency.



Monolithic Inserter
(Order by Request)



MIS Inserter



Modular Implant Inserter
Adjustable handle for straight, 45°, and 90° insertions



Corner Tamp




Articulating Tamp



Curette



Straight Trial



Rotating Cutter



Paddle Distractor

Discectomy Set 89-0031

Part #	Description	Qty
Top Tray		
89-0300	Anatomic Rasp	1
89-0301	Flat Rasp	1
89-0303	Curved Rasp	1
89-0600	Nerve Root Retractor, 70°	1
89-0601	Nerve Root Retractor, 90°	1
89-0602	Malleable Nerve Root Retractor	1
89-0650	2mm Pituitary Rongeur - Straight Biting	1
89-0651	2mm Pituitary Rongeur - Up Biting	1
89-0652	4mm Pituitary Rongeur - Straight Biting	1
89-0653	4mm Pituitary Rongeur - Up Biting	1

Middle Tray

89-0372	#0 Curette Straight- Serrated	1
89-0373	#1 Curette Straight- Serrated	1
89-0376	#0 Curette Up- Serrated	1
89-0377	#1 Curette Up- Serrated	1
89-0380	#0 Curette Down- Serrated	1
89-0381	#1 Curette Down- Serrated	1
89-0385	#1 Curette Right- Serrated	1
89-0389	#1 Curette Left- Serrated	1
89-0390	Small Ring Curette - 6W x 10L	1
89-0391	Large Ring Curette - 10W x 10L	1

Bottom Tray

89-0250	6mm Paddle Distractor	1
89-0251	7mm Paddle Distractor	1
89-0252	8mm Paddle Distractor	1
89-0253	9mm Paddle Distractor	1
89-0254	10mm Paddle Distractor	1
89-0255	11mm Paddle Distractor	1
89-0256	12mm Paddle Distractor	1
89-0257	13mm Paddle Distractor	1
89-0258	14mm Paddle Distractor	1
89-0259	15mm Paddle Distractor	1
89-0260	16mm Paddle Distractor	1


Part #	Description	Qty
89-0406	6mm Rotating Cutter	1
89-0407	7mm Rotating Cutter	1
89-0408	8mm Rotating Cutter	1
89-0409	9mm Rotating Cutter	1
89-0410	10mm Rotating Cutter	1
89-0411	11mm Rotating Cutter	1
89-0412	12mm Rotating Cutter	1
89-0413	13mm Rotating Cutter	1
89-0414	14mm Rotating Cutter	1
89-0415	15mm Rotating Cutter	1
89-0416	16mm Rotating Cutter	1

Discectomy Case


89-0092	Discectomy Case	
---------	-----------------	--

Order by Request - Discectomy Instruments


89-0351	#00 Curette Straight - Smooth
89-0353	#1 Curette Straight - Smooth
89-0355	#00 Curette Up - Smooth
89-0356	#0 Curette Up - Smooth
89-0357	#1 Curette Up - Smooth
89-0359	#00 Curette Down - Smooth
89-0360	#0 Curette Down - Smooth
89-0361	#1 Curette Down - Smooth
89-0364	#0 Curette Right - Smooth
89-0365	#1 Curette Right - Smooth
89-0368	#0 Curette Left - Smooth
89-0369	#1 Curette Left - Smooth
89-0342	#00 Curette 90° Down – Smooth
89-0343	#0 Curette 90° Down – Smooth
89-0344	#1 Curette 90° Down – Smooth
89-0346	#00 Curette 90° Down – Serrated
89-0347	#0 Curette 90° Down – Serrated
89-0348	#1 Curette 90° Down – Serrated



Top Level



Middle Level



Bottom Level

Instrument Set 89-0021


Part #	Description	Qty
Top Tray		
89-0100	Universal Handle	2
NA	9mm Straight Modular Implant Inserter (assembled)	
NA	11mm Straight Modular Implant Inserter (assembled)	
89-0113	Implant Inserter Body (C-Shaft) - (included within assembled Modular Implant Inserters)	2
89-0108	9mm Straight Implant Remover	1
89-0110	11mm Straight Implant Remover	1
89-0112	MIS Inserter	1
89-0114	9mm Straight Implant Shaft (1 is included within the assembled	2
89-0116	11mm Straight Implant Shaft	2
89-0117	9mm Curved Implant Shaft (1 is included within the assembled Modular Implant Inserter provided in the set)	2
89-0119	11mm Curved Implant Shaft	2
89-0120	Threaded Insert	6
89-0121	Impact Cap	1
89-0200	Slap Hammer	1
Bottom Tray		
89-0153	9mm Tamp	1
89-0155	11mm Tamp	1
89-0156	Corner (Y) Tamp	1
35-5001	Funnel	1
35-5002	Plunger	1
89-0500	Bone Packing Tool	Order By Request
89-0509	9mm Monolithic Straight Implant Inserter	Order By Request
89-0511	11mm Monolithic Straight Implant Inserter	Order By Request
89-0609	9mm Articulating Tamp	Order By Request
89-0611	11mm Articulating Tamp	Order By Request

Instruments Case


89-0094	Instruments Case	1
---------	------------------	---

Trial Set 89-0022 - Order By Request


Part #	Description	Qty
89-0093	Trial Case	1
89-0206	9W x 27L x 6H, 0° Straight Trial	1
89-0207	9W x 27L x 7H, 0° Straight Trial	1
89-0208	9W x 27L x 8H, 0° Straight Trial	1
89-0209	9W x 27L x 9H, 0° Straight Trial	1
89-0210	9W x 27L x 10H, 0° Straight Trial	1
89-0211	9W x 27L x 11H, 0° Straight Trial	1
89-0212	9W x 27L x 12H, 0° Straight Trial	1
89-0213	9W x 27L x 13H, 0° Straight Trial	1
89-0214	9W x 27L x 14H, 0° Straight Trial	1
89-0215	9W x 27L x 15H, 0° Straight Trial	1
89-0216	9W x 27L x 16H, 0° Straight Trial	1
89-0306	9W x 27L x 6H, 0° Curved Trial	1
89-0307	9W x 27L x 7H, 0° Curved Trial	1
89-0308	9W x 27L x 8H, 0° Curved Trial	1
89-0309	9W x 27L x 9H, 0° Curved Trial	1
89-0310	9W x 27L x 10H, 0° Curved Trial	1
89-0311	9W x 27L x 11H, 0° Curved Trial	1
89-0312	9W x 27L x 12H, 0° Curved Trial	1
89-0313	9W x 27L x 13H, 0° Curved Trial	1
89-0314	9W x 27L x 14H, 0° Curved Trial	1
89-0315	9W x 27L x 15H, 0° Curved Trial	1
89-0316	9W x 27L x 16H, 0° Curved Trial	1



One-Level



Top Level



Bottom Level

Straight Spacers

9mm x 23mm

89-1006SP	9W x 23L x 6, 0° Straight
89-1007SP	9W x 23L x 7, 0° Straight
89-1008SP	9W x 23L x 8, 0° Straight
89-1009SP	9W x 23L x 9, 0° Straight
89-1010SP	9W x 23L x 10, 0° Straight
89-1011SP	9W x 23L x 11, 0° Straight
89-1012SP	9W x 23L x 12, 0° Straight
89-1013SP	9W x 23L x 13, 0° Straight
89-1014SP	9W x 23L x 14, 0° Straight
89-1015SP	9W x 23L x 15, 0° Straight
89-1016SP	9W x 23L x 16, 0° Straight
89-7008SP	9W x 23L x 8, 8° Straight
89-7009SP	9W x 23L x 9, 8° Straight
89-7010SP	9W x 23L x 10, 8° Straight
89-7011SP	9W x 23L x 11, 8° Straight
89-7012SP	9W x 23L x 12, 8° Straight
89-7013SP	9W x 23L x 13, 8° Straight
89-7014SP	9W x 23L x 14, 8° Straight
89-7015SP	9W x 23L x 15, 8° Straight
89-7016SP	9W x 23L x 16, 8° Straight

9mm x 27mm

89-3006SP	9W x 27L x 6, 0° Straight
89-3007SP	9W x 27L x 7, 0° Straight
89-3008SP	9W x 27L x 8, 0° Straight
89-3009SP	9W x 27L x 9, 0° Straight
89-3010SP	9W x 27L x 10, 0° Straight
89-3011SP	9W x 27L x 11, 0° Straight
89-3012SP	9W x 27L x 12, 0° Straight
89-3013SP	9W x 27L x 13, 0° Straight
89-3014SP	9W x 27L x 14, 0° Straight
89-3015SP	9W x 27L x 15, 0° Straight
89-3016SP	9W x 27L x 16, 0° Straight
89-9009SP	9W x 27L x 9, 8° Straight
89-9010SP	9W x 27L x 10, 8° Straight
89-9011SP	9W x 27L x 11, 8° Straight
89-9012SP	9W x 27L x 12, 8° Straight
89-9013SP	9W x 27L x 13, 8° Straight
89-9014SP	9W x 27L x 14, 8° Straight
89-9015SP	9W x 27L x 15, 8° Straight
89-9016SP	9W x 27L x 16, 8° Straight

11mm x 23mm

89-2006SP	11W x 23L x 6, 0° Straight
89-2007SP	11W x 23L x 7, 0° Straight
89-2008SP	11W x 23L x 8, 0° Straight
89-2009SP	11W x 23L x 9, 0° Straight
89-2010SP	11W x 23L x 10, 0° Straight
89-2011SP	11W x 23L x 11, 0° Straight
89-2012SP	11W x 23L x 12, 0° Straight
89-2013SP	11W x 23L x 13, 0° Straight
89-2014SP	11W x 23L x 14, 0° Straight
89-2015SP	11W x 23L x 15, 0° Straight
89-2016SP	11W x 23L x 16, 0° Straight
89-8008SP	11W x 23L x 8, 8° Straight
89-8009SP	11W x 23L x 9, 8° Straight
89-8010SP	11W x 23L x 10, 8° Straight
89-8011SP	11W x 23L x 11, 8° Straight
89-8012SP	11W x 23L x 12, 8° Straight
89-8013SP	11W x 23L x 13, 8° Straight
89-8014SP	11W x 23L x 14, 8° Straight
89-8015SP	11W x 23L x 15, 8° Straight
89-8016SP	11W x 23L x 16, 8° Straight

11mm x 27mm


89-4006SP	11W x 27L x 6, 0° Straight
89-4007SP	11W x 27L x 7, 0° Straight
89-4008SP	11W x 27L x 8, 0° Straight
89-4009SP	11W x 27L x 9, 0° Straight
89-4010SP	11W x 27L x 10, 0° Straight
89-4011SP	11W x 27L x 11, 0° Straight
89-4012SP	11W x 27L x 12, 0° Straight
89-4013SP	11W x 27L x 13, 0° Straight
89-4014SP	11W x 27L x 14, 0° Straight
89-4015SP	11W x 27L x 15, 0° Straight
89-4016SP	11W x 27L x 16, 0° Straight
89-1209SP	11W x 27L x 9, 8° Straight
89-1210SP	11W x 27L x 10, 8° Straight
89-1211SP	11W x 27L x 11, 8° Straight
89-1212SP	11W x 27L x 12, 8° Straight
89-1213SP	11W x 27L x 13, 8° Straight
89-1214SP	11W x 27L x 14, 8° Straight
89-1215SP	11W x 27L x 15, 8° Straight
89-1216SP	11W x 27L x 16, 8° Straight

11mm x 31mm

89-6006SP	11W x 31L x 6, 0° Straight
89-6007SP	11W x 31L x 7, 0° Straight
89-6008SP	11W x 31L x 8, 0° Straight
89-6009SP	11W x 31L x 9, 0° Straight
89-6010SP	11W x 31L x 10, 0° Straight
89-6011SP	11W x 31L x 11, 0° Straight
89-6012SP	11W x 31L x 12, 0° Straight
89-6013SP	11W x 31L x 13, 0° Straight
89-6014SP	11W x 31L x 14, 0° Straight
89-6015SP	11W x 31L x 15, 0° Straight
89-6016SP	11W x 31L x 16, 0° Straight
89-1409SP	11W x 31L x 9, 8° Straight
89-1410SP	11W x 31L x 10, 8° Straight
89-1411SP	11W x 31L x 11, 8° Straight
89-1412SP	11W x 31L x 12, 8° Straight
89-1413SP	11W x 31L x 13, 8° Straight
89-1414SP	11W x 31L x 14, 8° Straight
89-1415SP	11W x 31L x 15, 8° Straight
89-1416SP	11W x 31L x 16, 8° Straight

11mm x 33mm

88-2006SP	11W x 33L x 6H, 0° Straight
88-2007SP	11W x 33L x 7H, 0° Straight
88-2008SP	11W x 33L x 8H, 0° Straight
88-2009SP	11W x 33L x 9H, 0° Straight
88-2010SP	11W x 33L x 10H, 0° Straight
88-2011SP	11W x 33L x 11H, 0° Straight
88-2012SP	11W x 33L x 12H, 0° Straight
88-2013SP	11W x 33L x 13H, 0° Straight
88-2014SP	11W x 33L x 14H, 0° Straight
88-2015SP	11W x 33L x 15H, 0° Straight
88-2016SP	11W x 33L x 16H, 0° Straight



9mm Straight Spacer 0°

Curved Spacers

9mm x 27mm

89-1706SP	9W x 27L x 6H, 0° Curved
89-1707SP	9W x 27L x 7H, 0° Curved
89-1708SP	9W x 27L x 8H, 0° Curved
89-1709SP	9W x 27L x 9H, 0° Curved
89-1710SP	9W x 27L x 10H, 0° Curved
89-1711SP	9W x 27L x 11H, 0° Curved
89-1712SP	9W x 27L x 12H, 0° Curved
89-1713SP	9W x 27L x 13H, 0° Curved
89-1714SP	9W x 27L x 14H, 0° Curved
89-1715SP	9W x 27L x 15H, 0° Curved
89-1716SP	9W x 27L x 16H, 0° Curved
89-4508SP	9W x 27L x 8H, 8° Curved
89-4509SP	9W x 27L x 9H, 8° Curved
89-4510SP	9W x 27L x 10H, 8° Curved
89-4511SP	9W x 27L x 11H, 8° Curved
89-4512SP	9W x 27L x 12H, 8° Curved
89-4513SP	9W x 27L x 13H, 8° Curved
89-4514SP	9W x 27L x 14H, 8° Curved
89-4515SP	9W x 27L x 15H, 8° Curved
89-4516SP	9W x 27L x 16H, 8° Curved

9mm X 31mm


89-1906SP	9W x 31L x 6H, 0° Curved
89-1907SP	9W x 31L x 7H, 0° Curved
89-1908SP	9W x 31L x 8H, 0° Curved
89-1909SP	9W x 31L x 9H, 0° Curved
89-1910SP	9W x 31L x 10H, 0° Curved
89-1911SP	9W x 31L x 11H, 0° Curved
89-1912SP	9W x 31L x 12H, 0° Curved
89-1913SP	9W x 31L x 13H, 0° Curved
89-1914SP	9W x 31L x 14H, 0° Curved
89-1915SP	9W x 31L x 15H, 0° Curved
89-1916SP	9W x 31L x 16H, 0° Curved
89-4708SP	9W x 31L x 8H, 8° Curved
89-4709SP	9W x 31L x 9H, 8° Curved
89-4710SP	9W x 31L x 10H, 8° Curved
89-4711SP	9W x 31L x 11H, 8° Curved
89-4712SP	9W x 31L x 12H, 8° Curved
89-4713SP	9W x 31L x 13H, 8° Curved
89-4714SP	9W x 31L x 14H, 8° Curved
89-4715SP	9W x 31L x 15H, 8° Curved
89-4716SP	9W x 31L x 16H, 8° Curved

11mm x 27mm

89-1806SP	11W x 27L x 6H, 0° Curved
89-1807SP	11W x 27L x 7H, 0° Curved
89-1808SP	11W x 27L x 8H, 0° Curved
89-1809SP	11W x 27L x 9H, 0° Curved
89-1810SP	11W x 27L x 10H, 0° Curved
89-1811SP	11W x 27L x 11H, 0° Curved
89-1812SP	11W x 27L x 12H, 0° Curved
89-1813SP	11W x 27L x 13H, 0° Curved
89-1814SP	11W x 27L x 14H, 0° Curved
89-1815SP	11W x 27L x 15H, 0° Curved
89-1816SP	11W x 27L x 16H, 0° Curved
89-4608SP	11W x 27L x 8H, 8° Curved
89-4609SP	11W x 27L x 9H, 8° Curved
89-4610SP	11W x 27L x 10H, 8° Curved
89-4611SP	11W x 27L x 11H, 8° Curved
89-4612SP	11W x 27L x 12H, 8° Curved
89-4613SP	11W x 27L x 13H, 8° Curved
89-4614SP	11W x 27L x 14H, 8° Curved
89-4615SP	11W x 27L x 15H, 8° Curved
89-4616SP	11W x 27L x 16H, 8° Curved

11mm X 31mm

89-4206SP	11W x 31L x 6H, 0° Curved
89-4207SP	11W x 31L x 7H, 0° Curved
89-4208SP	11W x 31L x 8H, 0° Curved
89-4209SP	11W x 31L x 9H, 0° Curved
89-4210SP	11W x 31L x 10H, 0° Curved
89-4211SP	11W x 31L x 11H, 0° Curved
89-4212SP	11W x 31L x 12H, 0° Curved
89-4213SP	11W x 31L x 13H, 0° Curved
89-4214SP	11W x 31L x 14H, 0° Curved
89-4215SP	11W x 31L x 15H, 0° Curved
89-4216SP	11W x 31L x 16H, 0° Curved
89-4808SP	11W x 31L x 8H, 8° Curved
89-4809SP	11W x 31L x 9H, 8° Curved
89-4810SP	11W x 31L x 10H, 8° Curved
89-4811SP	11W x 31L x 11H, 8° Curved
89-4812SP	11W x 31L x 12H, 8° Curved
89-4813SP	11W x 31L x 13H, 8° Curved
89-4814SP	11W x 31L x 14H, 8° Curved
89-4815SP	11W x 31L x 15H, 8° Curved
89-4816SP	11W x 31L x 16H, 8° Curved



9mm Curved Spacer 0°

