

In a world of trauma  
a **GALAXY**<sup>TM</sup> of solutions




## HUMERUS

N.	Code	Description
3	99-93100	Threaded Wire Sterile (pack of 2)
3	93620	Wire Locking Clamp
4	93310	Small Clamp
1	936080	Rod Ø 6mm L 80mm
1	936140	Rod Ø 6mm L 140mm




N.	Code
6	93110
2	939100
1	939150
4	911530




## WRIST

N.	Code	Description
5	93310	Small Clamp
1	936140	Rod Ø 6mm L 140mm
1	936080	Rod Ø 6mm L 80mm
4	947320	Self Drilling XCaliber Cylindrical Screw Shaft Ø 4mm Thread 3mm L 100/20 QC




N.	Code
6	93310
1	936100
1	936080
1	936060
4	947320



## PELVIS

N.	Code	Description
8	93010	Large Clamp
2	932200	Rod Ø 12mm L 200mm
1	932300	Rod Ø 12mm L 300mm
1	932400	Rod Ø 12mm L 400mm
4	912650	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 260/50



## FEMUR

N.	Code
6	93010
1	932100
2	932200
4	912640

# Simple components,

## Description

Medium Clamp  
Rod Ø 9mm L 100mm  
Rod Ø 9mm L 150mm  
Self-drilling XCaliber Screws  
Shaft Ø 6mm  
Thread Ø 6.0-5.6mm 150/30



## ELBOW

N.	Code	Description
1	93020	Multiscrew Clamp
1	932200	Rod Ø 12mm L 200mm
1	93410	Elbow Hinge
1	939150	Rod Ø 9mm L 150mm
2	93110	Medium Clamp
2	911530	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 150/30
2	10137	Cortical Screws Shaft Ø 6mm Thread Ø 4.5-3.5mm 120/20

## Description

Small Clamp  
Rod Ø 6mm L 100mm  
Rod Ø 6mm L 80mm  
Rod Ø 6mm L 60mm  
Self Drilling XCaliber Cylindrical  
Screw Shaft Ø 4mm  
Thread 3mm L 100/20 QC



N.	Code	Description
1	93350	Wrist Module
2	93330	Small Multiscrew Clamp-Short
4	947320	Self Drilling XCaliber Cylindrical Screw Shaft Ø 4mm Thread 3mm L 100/20 QC


## Description

Large Clamp  
Rod Ø 12mm L 100mm  
Rod Ø 12mm L 200mm  
XCaliber Screws  
260/40 6.0-5.6mm



## KNEE

N.	Code	Description
8	93010	Large Clamp
2	932350	Rod Ø 12mm L 350mm
1	932250	Rod Ø 12mm L 250mm
4	911540	XCaliber Screw 150/40 6-5.6mm
2	911550	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 150/50



## FOOT & ANKLE

N.	Code	Description
8	93010	Large Clamp
1	932200	Rod Ø 12mm L 200mm
1	932100	Rod Ø 12mm L 100mm
1	932250	Rod Ø 12mm L 250mm
1	932400	Rod Ø 12mm L 400mm
1	911540	XCaliber Screw 150/40 6-5.6mm
2	911530	XCaliber Screw 150/30 6-5.6mm


THE SYSTEM IS

# Infinite solutions



FOREARM

N.	Code	Description
5	93110	Medium Clamp
1	939100	Rod Ø 9mm L 100mm
1	939150	Rod Ø 9mm L 150mm
4	10137	Cortical Screws Shaft Ø 6mm Thread Ø 4.5-3.5mm 120/20



N.	Code	Description
1	936010	6mm L Rod
1	93330	Small Multiscrew Clamp-Short
1	93320	Small Multiscrew Clamp-Long
4	947320	Self Drilling XCaliber Cylindrical Screw Shaft Ø 4mm Thread 3mm L 100/20 QC




PELVIS



TIBIA

N.	Code	Description
6	93010	Large Clamp
1	932100	Rod Ø 12mm L 100mm
1	932150	Rod Ø 12mm L 150mm
1	932200	Rod Ø 12mm L 200mm
4	911540	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 150/40




N.	Code	Description
4	93010	Large Clamp
2	99-93030	Large-Medium Transition Clamp (Sterile)
1	932350	Rod Ø 12mm L 350mm
1	932250	Rod Ø 12mm L 250mm
1	932200	Rod Ø 12mm L 200mm
1	92080	Transfix Pin 80mm, Ø 4mm
2	911540	XCaliber Screw 150/40 6-5.6mm

Possibility to use six Large Clamps  
and Ø 6mm Transfix Pin (93080)



N.	Code	Description
6	93010	Large Clamp
2	932250	Rod Ø 12mm L 250mm
1	932350	Rod Ø 12mm L 350mm
4	911550	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 150/50



N.	Code	Description
9	93010	Large Clamp
3	932300	Rod Ø 12mm L 300mm
2	932250	Rod Ø 12mm L 250mm
4	911550	Self-drilling XCaliber Screws Shaft Ø 6mm Thread Ø 6.0-5.6mm 150/50



N.	Code	Description
1	99-93794	Galaxy UNYCO Diaphyseal Tibia Box
		or
1	99-93574	Galaxy UNYCO Mini Kit Tibia Sterile
2	99-932350	Rod D12mm L 350mm Sterile
1	99-93509	Galaxy UNYCO Mini Kit Instruments Sterile

# Components

The Galaxy Fixation components are MR conditional when marked with "MR".


## Clamps



## Multiscrew Clamps



## Modules



# Components

The Galaxy Fixation components are MR conditional when marked with "△/MR".

## Rods Diameter 12mm



Code	Description
932100	Rod 100mm long
932150	Rod 150mm long
932200	Rod 200mm long
932250	Rod 250mm long
932300	Rod 300mm long
932350	Rod 350mm long
932400	Rod 400mm long
99-932450	Rod 450mm long, sterile*
99-932500	Rod 500mm long, sterile*
99-932550	Rod 550mm long, sterile*
99-932600	Rod 600mm long, sterile*
99-932650	Rod 650mm long, sterile*

\* Special order only

## Rods Diameter 9mm



Code	Description
939100	Rod 100mm long
939150	Rod 150mm long
939200	Rod 200mm long
939250	Rod 250mm long
939300	Rod 300mm long

## Rods Diameter 6mm




Code	Description
936060	Rod 60mm long
936080	Rod 80mm long
936100	Rod 100mm long
936120	Rod 120mm long
936140	Rod 140mm long
936150	Rod 150mm long
936160	Rod 160mm long
936180	Rod 180mm long
936200	Rod 200mm long

## L Rod Diameter 6mm



Code	Description
936010	6mm L Rod

## Semi-Circular Rods Diam. 12mm (Aluminium)



Code	Description
932010	Semi-Circular Rod small 180mm long
932020	Semi-Circular Rod medium 215mm long
932030	Semi-Circular Rod large 250mm long

## Semi-Circular Rods Diam. 9mm

939010	Semi-Circular Rod small 115mm long
939020	Semi-Circular Rod medium 140mm long
939030	Semi-circular Rod large 165mm long

**Note:** All components of Galaxy Fixation System frames must be identified as MR Conditional prior to being placed in or near an MR Environment.

The Threaded Wires (93100), the Wire Locking Clamps (93620), and the L Rod (936010) and Semi-Circular Rods (939010, 939020, 939030) are not MR Conditional. Any construct/frame that is using Threaded Wires, the Wire Locking Clamps, the L-Rod and Semi-Circular Rods must therefore be considered as MR Unsafe.

THE SYSTEM IS MRI CONDITIONAL at 1.5 and 3 Tesla  
More details are available on relevant Instructions  
For Use - PQGAL



### Orthofix ABS – A touch of Orthofix Knowledge

We are now covering all the mobile devices:  
iOS version (App Store) for iPad  
Android version (Google Play Store) for Android tablet  
Website: <http://abs.orthofix.it/home.asp> for smartphone and PC



Try the app and then insert your review!  
This will help us to always innovate.


# Get it all sterile!

Besides the sterile pre-packaged kits, Galaxy Fixation System offers all **clamps and rods single-packed and sterile**, except for the 6mm L-rod cod. 936010. They can be ordered using the item code number preceded by 99- (e.g. 99-939300).


**99-93505** Galaxy Shoulder Sterile Kit



**99-93504** Galaxy Elbow Sterile Kit




**99-93501** Pelvis Sterile Kit




**99-93498** Galaxy Small Z-Configuration Sterile Kit



**99-93502** Lower Limb Diaphyseal Sterile Kit




**99-93499** Ankle Sterile kit - Transfix Pin  $\varnothing$ 4mm




**99-93503** Ankle Sterile Kit




**99-93520** Galaxy Medium - Paediatric Sterile Kit -  $\varnothing$ 5mm screw thread



**99-93601** Galaxy Wrist Sterile Kit



**99-93521** Galaxy Medium - Paediatric Sterile Kit -  $\varnothing$ 4mm screw thread




Manufactured by:  
ORTHOFIX Srl  
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy  
Telephone +39 045 6719000, Fax +39 045 6719380

[www.orthofix.com](http://www.orthofix.com)



GF-1202-PL-E0 M 01/19

Distributed by:



  
**ORTHOFIX®**