



CHIMAERA - Hip Fracture System, Instrument Tray

Top Tray




Mid Tray



Instrument Tray			
Code		Description	Q.ty
193990C		CHIMAERA - Hip Fracture System, Instrument tray, Complete	
193990		CHIMAERA - Hip Fracture System, Instrument tray, Empty	
TOP TRAY			
193115	1	Impaction Rod	1
193110	2	Locking Bolt	1
193270	3	Entry Reamer (max. Ø 15.5mm)	1
193260	4	Cannulated Awl	1
193230	5	Entry Reamer Guide	1
193325	6	6mm Hex Screwdriver	1
193100	7	Targeting Handle	1
91017	8	Universal Allen Wrench	1
193271	9	6mm Hex Extension	1
177380	10	Slotted Mallet	1
17955	11	Universal Chuck With T-Handle	1
193275	12	Sizing Guide	1
MID TRAY			
193320	1	Distal Screwdriver	1
193319	2	Distal Screwdriver Retention Rod	1
193211	3	Distal Tissue Guide	1
193212	4	Distal Trocar	1
193231	5	Lag Screw Tissue Guide	1
193232	6	Lag Screw Trocar	1
193274	7	Screw Sizing Gauge	1
193283	8	Lag Screwdriver	1
193282	9	Lag Screwdriver Retention Rod	1
193286	10	4,2mm Distal Drill Bit - Long	1
193281	11	Lag Screw Compression Rod	1
193970	12	Lag Screw Gradual Reamer	1
193973	13	4mm Cortical Drill	1
193321	14	Distal Screwdriver - Short	1
193318	15	Distal Screwdriver Retention Rod - Short	1
193285	16	4,2mm Distal Drill Bit - Short	2
193213	17	Distal Tissue Guide - Short	1
193277	18	Distal Screw Sizing Gauge	1
193284	19	Wire Pusher	1

CHIMAERA - Hip Fracture System, Instrument Tray

Base Tray




BASE TRAY

193222	1	Supplementary Screw Trocar	1
193221	2	Supplementary Tissue Guide	1
193292	3	Supplementary Screwdriver Retention Rod	1
193293	4	Supplementary Screwdriver	1
193971	5	Supplementary Gradual Drill Bit	1
193279	6	Hj Quick Connect T Handle	1
193276	7	Long Nail Ruler	1
173276	8	Ruler Support	1
193948	9	Multi Hole Wire Guide	1

Out of tray

Code	Description
99-193287	Threaded 3.2mm guide wire - 400mm - Sterile
172001	Flexible Reaming System Instruments Box, Complete
99-173281	Guide Wire with olive Ø 3x980mm
193265	Implant Template

CHIMAERA - Hip Fracture System, Extraction Instrument Tray



Extraction Instrument Tray		
Code	Description	Q.ty
193991C	CHIMAERA - Hip Fracture System, Extraction Instrument Tray, Complete	
193991	CHIMAERA - Hip Fracture System, Extraction Instrument Tray, Empty	
SMN173370	1 Slap Hammer	1
177380	2 Slotted Mallet	1
193325	3 6mm Hex Screwdriver	1
193271	4 6mm Hex Extension	1
193283	5 Lag Screwdriver	1
193280	6 Lag Screwdriver Retention Rod - Short	1
193293	7 Supplementary Screwdriver	1
193290	8 Supplementary Screwdriver Retention Rod - Short	1
193321	9 Distal Screwdriver - Short	1
193318	10 Distal Screwdriver Retention Rod - Short	1
17955	11 Universal Chuck With T-Handle	1
17978	12 Plier Cutter	1
177395	13 Needle Nose Vise Grips	1
193336	14 M8 Conical Mallet Adaptor	1
193337	15 Screw Extractor Size 4 - 9mm	1
91017	16 Universal Allen Wrench	1