


ContoursLPS

LAPIDUS PLATING SYSTEM



CONTOURED SPECIFICALLY FOR THE ANATOMY



Precisely contoured for the first TMT



Plate Thickness


Contours LPS: 1.78mm (avg)

Comparative competitive thickness = 1.55mm



Three points of fixation proximal and distal to the arthrodesis site

Locking screw angles are designed to capture the "sweet" spots of the bone – allowing for optimal purchase



Patent pending trapezoidal screw configuration for fusion stability

2mm of uniform, mechanical compression





Pre-Op



Post-Op



Pre-Op



Post-Op



ORTHOFIX: REDEFINING YOUR OPTIONS


Trinity® Evolution™
An Allograft with
Viable Cells



Centronail
Ankle Compression
Nailing System



Collage™
Osteoconductive
Scaffold



TrueLok
Ring Fixation
System

Physio-Stim®
Bone Growth
Stimulation



MiniRail™
Fixators



The Guided Growth™ System
eight-Plate®



Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000, Fax +39 045 6719380

www.orthofix.com

LP-1307-PL-E0 B 02/16


ORTHOFIX®