

1 HUM	2 HUM - TI	3 ELB	4 FA - WR - PAE																																																																																																												
99-93505 Galaxy Shoulder Sterile Kit: <table border="1"> <thead> <tr> <th>Code</th> <th>Description</th> </tr> </thead> <tbody> <tr><td>4x93310</td><td>Small Clamp</td></tr> <tr><td>3x93620</td><td>Wire Locking Clamp</td></tr> <tr><td>1x936080</td><td>Rod d 6mm L 80mm</td></tr> <tr><td>1x936100</td><td>Rod d 6mm L 100mm</td></tr> <tr><td>1x936120</td><td>Rod d 6mm L 120mm</td></tr> <tr><td>1x936140</td><td>Rod d 6mm L 140mm</td></tr> <tr><td>1x30017</td><td>Allen Wrench 5mm</td></tr> <tr><td>1x19970</td><td>Wire Guide</td></tr> <tr><td>1x91150</td><td>Universal "T" Wrench</td></tr> <tr><td>1x81031</td><td>Open end wrench</td></tr> <tr><td>6x93100</td><td>300mm Threaded wire</td></tr> </tbody> </table>	Code	Description	4x93310	Small Clamp	3x93620	Wire Locking Clamp	1x936080	Rod d 6mm L 80mm	1x936100	Rod d 6mm L 100mm	1x936120	Rod d 6mm L 120mm	1x936140	Rod d 6mm L 140mm	1x30017	Allen Wrench 5mm	1x19970	Wire Guide	1x91150	Universal "T" Wrench	1x81031	Open end wrench	6x93100	300mm Threaded wire	99-91600UE XCaliber Metadiaphyseal Complete Kit, sterile: <table border="1"> <thead> <tr> <th>Code</th> <th>Description</th> </tr> </thead> <tbody> <tr><td>1x91150</td><td>Universal T Wrench</td></tr> <tr><td>1x91017</td><td>Universal Allen Wrench</td></tr> <tr><td>1x11106</td><td>Drill Guide 3.2x40mm</td></tr> <tr><td>1x11104</td><td>Drill Guide 4.8x40mm</td></tr> <tr><td>1x91000</td><td>XCaliber Meta-Diaphyseal Fixator</td></tr> <tr><td>3x11102</td><td>Screw Guide L 60mm</td></tr> <tr><td>1x91015</td><td>Dynamic Compression/Distraction Unit</td></tr> <tr><td>1x99-611560</td><td>XCaliber Osteotite Screw 150/60</td></tr> <tr><td>1x99-611550</td><td>XCaliber Osteotite Screw 150/50</td></tr> <tr><td>1x99-611540</td><td>XCaliber Osteotite Screw 150/40</td></tr> <tr><td>3x99-611530</td><td>XCaliber Osteotite Screws 150/30</td></tr> <tr><td>1x1101101</td><td>Cannulated Drill Bit 3.2x200mm</td></tr> <tr><td>2x110114</td><td>X-Wires without olive 1.5x250mm</td></tr> <tr><td>1x1100101</td><td>Drill Bit 4.8x180mm</td></tr> </tbody> </table>	Code	Description	1x91150	Universal T Wrench	1x91017	Universal Allen Wrench	1x11106	Drill Guide 3.2x40mm	1x11104	Drill Guide 4.8x40mm	1x91000	XCaliber Meta-Diaphyseal Fixator	3x11102	Screw Guide L 60mm	1x91015	Dynamic Compression/Distraction Unit	1x99-611560	XCaliber Osteotite Screw 150/60	1x99-611550	XCaliber Osteotite Screw 150/50	1x99-611540	XCaliber Osteotite Screw 150/40	3x99-611530	XCaliber Osteotite Screws 150/30	1x1101101	Cannulated Drill Bit 3.2x200mm	2x110114	X-Wires without olive 1.5x250mm	1x1100101	Drill Bit 4.8x180mm	99-93504 Galaxy Elbow Sterile Kit: <table border="1"> <thead> <tr> <th>Code</th> <th>Description</th> </tr> </thead> <tbody> <tr><td>1x93020</td><td>Multiscrew Clamp</td></tr> <tr><td>2x93110</td><td>Medium Clamp</td></tr> <tr><td>1x93410</td><td>Elbow Hinge</td></tr> <tr><td>1x932200</td><td>Rod Ø 12mm L 200mm</td></tr> <tr><td>1x939150</td><td>Rod Ø 9mm L 150mm</td></tr> <tr><td>1x30017</td><td>Allen Wrench 5mm</td></tr> <tr><td>1x91150</td><td>Universal "T" Wrench</td></tr> <tr><td>1x1-1100201</td><td>Drill Bit Ø 4.8x240mm</td></tr> <tr><td>1x11138</td><td>Drill guide Ø 4.8mm</td></tr> <tr><td>2x11137</td><td>Screw Guide 80mm</td></tr> <tr><td>1x19950</td><td>Drill Guide Ø 3.2mm L 100mm</td></tr> <tr><td>1x1-1100301</td><td>Drill Bit Ø 3.2x200mm QC</td></tr> <tr><td>2x911530</td><td>Xcaliber Screws 150/30</td></tr> <tr><td>2x10137</td><td>Cortical Screws 120/20 Ø 4,5/3,5 Shaft 6mm</td></tr> <tr><td>1x11146</td><td>X-Wire without olive Ø 2mm Length 150mm</td></tr> </tbody> </table>	Code	Description	1x93020	Multiscrew Clamp	2x93110	Medium Clamp	1x93410	Elbow Hinge	1x932200	Rod Ø 12mm L 200mm	1x939150	Rod Ø 9mm L 150mm	1x30017	Allen Wrench 5mm	1x91150	Universal "T" Wrench	1x1-1100201	Drill Bit Ø 4.8x240mm	1x11138	Drill guide Ø 4.8mm	2x11137	Screw Guide 80mm	1x19950	Drill Guide Ø 3.2mm L 100mm	1x1-1100301	Drill Bit Ø 3.2x200mm QC	2x911530	Xcaliber Screws 150/30	2x10137	Cortical Screws 120/20 Ø 4,5/3,5 Shaft 6mm	1x11146	X-Wire without olive Ø 2mm Length 150mm	99-93498 Galaxy Small Z-Configuration Sterile Kit: <table border="1"> <thead> <tr> <th>Code</th> <th>Description</th> </tr> </thead> <tbody> <tr><td>6x93310</td><td>Small clamp</td></tr> <tr><td>1x936120</td><td>Rod Ø 6mm L 120mm</td></tr> <tr><td>1x936100</td><td>Rod Ø 6mm L 100mm</td></tr> <tr><td>1x936080</td><td>Rod Ø 6mm L 80mm</td></tr> <tr><td>2x19995</td><td>Sleeve 4.5mm</td></tr> <tr><td>2x13530</td><td>Drill guide d 2.7mm</td></tr> <tr><td>1x91017</td><td>Allen Wrench</td></tr> <tr><td>4x947320</td><td>Self Drilling XCaliber Cylindrical Screw Shaft Ø 4mm Thread 3mm L 100/20 QC</td></tr> <tr><td>2x1-1355001</td><td>Drill Bit d 2,7mm</td></tr> <tr><td>1x93160</td><td>"T" Wrench AO QC</td></tr> </tbody> </table>	Code	Description	6x93310	Small clamp	1x936120	Rod Ø 6mm L 120mm	1x936100	Rod Ø 6mm L 100mm	1x936080	Rod Ø 6mm L 80mm	2x19995	Sleeve 4.5mm	2x13530	Drill guide d 2.7mm	1x91017	Allen Wrench	4x947320	Self Drilling XCaliber Cylindrical Screw Shaft Ø 4mm Thread 3mm L 100/20 QC	2x1-1355001	Drill Bit d 2,7mm	1x93160	"T" Wrench AO QC
Code	Description																																																																																																														
4x93310	Small Clamp																																																																																																														
3x93620	Wire Locking Clamp																																																																																																														
1x936080	Rod d 6mm L 80mm																																																																																																														
1x936100	Rod d 6mm L 100mm																																																																																																														
1x936120	Rod d 6mm L 120mm																																																																																																														
1x936140	Rod d 6mm L 140mm																																																																																																														
1x30017	Allen Wrench 5mm																																																																																																														
1x19970	Wire Guide																																																																																																														
1x91150	Universal "T" Wrench																																																																																																														
1x81031	Open end wrench																																																																																																														
6x93100	300mm Threaded wire																																																																																																														
Code	Description																																																																																																														
1x91150	Universal T Wrench																																																																																																														
1x91017	Universal Allen Wrench																																																																																																														
1x11106	Drill Guide 3.2x40mm																																																																																																														
1x11104	Drill Guide 4.8x40mm																																																																																																														
1x91000	XCaliber Meta-Diaphyseal Fixator																																																																																																														
3x11102	Screw Guide L 60mm																																																																																																														
1x91015	Dynamic Compression/Distraction Unit																																																																																																														
1x99-611560	XCaliber Osteotite Screw 150/60																																																																																																														
1x99-611550	XCaliber Osteotite Screw 150/50																																																																																																														
1x99-611540	XCaliber Osteotite Screw 150/40																																																																																																														
3x99-611530	XCaliber Osteotite Screws 150/30																																																																																																														
1x1101101	Cannulated Drill Bit 3.2x200mm																																																																																																														
2x110114	X-Wires without olive 1.5x250mm																																																																																																														
1x1100101	Drill Bit 4.8x180mm																																																																																																														
Code	Description																																																																																																														
1x93020	Multiscrew Clamp																																																																																																														
2x93110	Medium Clamp																																																																																																														
1x93410	Elbow Hinge																																																																																																														
1x932200	Rod Ø 12mm L 200mm																																																																																																														
1x939150	Rod Ø 9mm L 150mm																																																																																																														
1x30017	Allen Wrench 5mm																																																																																																														
1x91150	Universal "T" Wrench																																																																																																														
1x1-1100201	Drill Bit Ø 4.8x240mm																																																																																																														
1x11138	Drill guide Ø 4.8mm																																																																																																														
2x11137	Screw Guide 80mm																																																																																																														
1x19950	Drill Guide Ø 3.2mm L 100mm																																																																																																														
1x1-1100301	Drill Bit Ø 3.2x200mm QC																																																																																																														
2x911530	Xcaliber Screws 150/30																																																																																																														
2x10137	Cortical Screws 120/20 Ø 4,5/3,5 Shaft 6mm																																																																																																														
1x11146	X-Wire without olive Ø 2mm Length 150mm																																																																																																														
Code	Description																																																																																																														
6x93310	Small clamp																																																																																																														
1x936120	Rod Ø 6mm L 120mm																																																																																																														
1x936100	Rod Ø 6mm L 100mm																																																																																																														
1x936080	Rod Ø 6mm L 80mm																																																																																																														
2x19995	Sleeve 4.5mm																																																																																																														
2x13530	Drill guide d 2.7mm																																																																																																														
1x91017	Allen Wrench																																																																																																														
4x947320	Self Drilling XCaliber Cylindrical Screw Shaft Ø 4mm Thread 3mm L 100/20 QC																																																																																																														
2x1-1355001	Drill Bit d 2,7mm																																																																																																														
1x93160	"T" Wrench AO QC																																																																																																														

5 WR	6 WR	7 PEL	8 FE - KN - TI																																																																																		
99-93601 Galaxy Wrist Sterile Kit: <table border="1"> <thead> <tr> <th>Code</th> <th>Description</th> </tr> </thead> <tbody> <tr><td>2x93330</td><td>Small multiscrew clamp-SHORT</td></tr> <tr><td>1x93350</td><td>Wrist Module</td></tr> <tr><td>2x19995</td><td>Screw Guide</td></tr> <tr><td>2x13530</td><td>Drill guide Ø 2.7mm</td></tr> <tr><td>1x91017</td><td>Allen Wrench</td></tr> <tr><td>4x947320</td><td>Self Drilling XCaliber Cylindrical Screw Shaft Ø 4mm Thread 3mm L 100/20 QC</td></tr> <tr><td>2x1-1355001</td><td>Drill Bit Ø 2,7mm</td></tr> <tr><td>1x93160</td><td>Screw T Wrench QC</td></tr> <tr><td>3x13715</td><td>Kwire Ø 1,5x150mm</td></tr> </tbody> </table>	Code	Description	2x93330	Small multiscrew clamp-SHORT	1x93350	Wrist Module	2x19995	Screw Guide	2x13530	Drill guide Ø 2.7mm	1x91017	Allen Wrench	4x947320	Self Drilling XCaliber Cylindrical Screw Shaft Ø 4mm Thread 3mm L 100/20 QC	2x1-1355001	Drill Bit Ø 2,7mm	1x93160	Screw T Wrench QC	3x13715	Kwire Ø 1,5x150mm	99-36501 RadioLucent Wrist Fixator Kit, sterile: <table border="1"> <thead> <tr> <th>Q.ty</th> <th>Description</th> </tr> </thead> <tbody> <tr><td>1x</td><td>RadioLucent Wrist Fixator</td></tr> <tr><td>1x</td><td>Compression-Distraction Unit</td></tr> <tr><td>1x</td><td>Allen Wrench 4mm</td></tr> <tr><td>2x</td><td>Drill Guides</td></tr> <tr><td>1x</td><td>Drill Bit Ø 2.7mm</td></tr> <tr><td>1x</td><td>Screw Template</td></tr> <tr><td>1x</td><td>T-Wrench</td></tr> <tr><td>2x</td><td>Self-drilling Screws, shaft Ø 4mm, thread Ø 3.3-3mm, total length 70mm, thread length 20mm</td></tr> <tr><td>2x</td><td>Self-drilling Screws, shaft Ø 4mm, thread Ø 3.3-3mm, total length 80mm, thread length 25mm</td></tr> <tr><td>4x</td><td>Bone Screw Covers</td></tr> </tbody> </table>	Q.ty	Description	1x	RadioLucent Wrist Fixator	1x	Compression-Distraction Unit	1x	Allen Wrench 4mm	2x	Drill Guides	1x	Drill Bit Ø 2.7mm	1x	Screw Template	1x	T-Wrench	2x	Self-drilling Screws, shaft Ø 4mm, thread Ø 3.3-3mm, total length 70mm, thread length 20mm	2x	Self-drilling Screws, shaft Ø 4mm, thread Ø 3.3-3mm, total length 80mm, thread length 25mm	4x	Bone Screw Covers	99-93501 Pelvis Sterile Kit: <table border="1"> <thead> <tr> <th>Code</th> <th>Description</th> </tr> </thead> <tbody> <tr><td>8x93010</td><td>Large Clamp</td></tr> <tr><td>1x932350</td><td>Rod Ø 12mm L 350mm</td></tr> <tr><td>1x932300</td><td>Rod Ø 12mm L 300mm</td></tr> <tr><td>2x932200</td><td>Rod Ø 12mm L 200mm</td></tr> <tr><td>4x912640</td><td>Self-drilling XCaliber Screws, L 260mm, thread length 40mm</td></tr> <tr><td>1x11138</td><td>Drill Guide Ø 4.8mm</td></tr> <tr><td>1x11137</td><td>Screw Guide 80mm</td></tr> <tr><td>1x1-1100101</td><td>Drill Bit Ø 4.8x180mm Coated - Quick Connect</td></tr> <tr><td>1x91150</td><td>Bone Screw T Wrench</td></tr> </tbody> </table>	Code	Description	8x93010	Large Clamp	1x932350	Rod Ø 12mm L 350mm	1x932300	Rod Ø 12mm L 300mm	2x932200	Rod Ø 12mm L 200mm	4x912640	Self-drilling XCaliber Screws, L 260mm, thread length 40mm	1x11138	Drill Guide Ø 4.8mm	1x11137	Screw Guide 80mm	1x1-1100101	Drill Bit Ø 4.8x180mm Coated - Quick Connect	1x91150	Bone Screw T Wrench	99-93502 Lower Limb Diaphyseal Sterile Kit: <table border="1"> <thead> <tr> <th>Code</th> <th>Description</th> </tr> </thead> <tbody> <tr><td>6x93010</td><td>Large Clamp</td></tr> <tr><td>1x932300</td><td>Rod Ø 12mm L 300mm</td></tr> <tr><td>2x932150</td><td>Rod Ø 12mm L 150mm</td></tr> <tr><td>2x912650</td><td>Self-drilling XCaliber Screws, L 260mm, thread length 50mm</td></tr> <tr><td>2x911550</td><td>Self-drilling XCaliber Screws, L 150mm, thread length 50mm</td></tr> <tr><td>1x11138</td><td>Drill Guide Ø 4.8mm</td></tr> <tr><td>1x11137</td><td>Screw Guide 80mm</td></tr> <tr><td>1x1-1100101</td><td>Drill Bit Ø 4.8x180mm Coated - Quick Connect</td></tr> <tr><td>1x91150</td><td>Bone Screw T Wrench</td></tr> </tbody> </table>	Code	Description	6x93010	Large Clamp	1x932300	Rod Ø 12mm L 300mm	2x932150	Rod Ø 12mm L 150mm	2x912650	Self-drilling XCaliber Screws, L 260mm, thread length 50mm	2x911550	Self-drilling XCaliber Screws, L 150mm, thread length 50mm	1x11138	Drill Guide Ø 4.8mm	1x11137	Screw Guide 80mm	1x1-1100101	Drill Bit Ø 4.8x180mm Coated - Quick Connect	1x91150	Bone Screw T Wrench
Code	Description																																																																																				
2x93330	Small multiscrew clamp-SHORT																																																																																				
1x93350	Wrist Module																																																																																				
2x19995	Screw Guide																																																																																				
2x13530	Drill guide Ø 2.7mm																																																																																				
1x91017	Allen Wrench																																																																																				
4x947320	Self Drilling XCaliber Cylindrical Screw Shaft Ø 4mm Thread 3mm L 100/20 QC																																																																																				
2x1-1355001	Drill Bit Ø 2,7mm																																																																																				
1x93160	Screw T Wrench QC																																																																																				
3x13715	Kwire Ø 1,5x150mm																																																																																				
Q.ty	Description																																																																																				
1x	RadioLucent Wrist Fixator																																																																																				
1x	Compression-Distraction Unit																																																																																				
1x	Allen Wrench 4mm																																																																																				
2x	Drill Guides																																																																																				
1x	Drill Bit Ø 2.7mm																																																																																				
1x	Screw Template																																																																																				
1x	T-Wrench																																																																																				
2x	Self-drilling Screws, shaft Ø 4mm, thread Ø 3.3-3mm, total length 70mm, thread length 20mm																																																																																				
2x	Self-drilling Screws, shaft Ø 4mm, thread Ø 3.3-3mm, total length 80mm, thread length 25mm																																																																																				
4x	Bone Screw Covers																																																																																				
Code	Description																																																																																				
8x93010	Large Clamp																																																																																				
1x932350	Rod Ø 12mm L 350mm																																																																																				
1x932300	Rod Ø 12mm L 300mm																																																																																				
2x932200	Rod Ø 12mm L 200mm																																																																																				
4x912640	Self-drilling XCaliber Screws, L 260mm, thread length 40mm																																																																																				
1x11138	Drill Guide Ø 4.8mm																																																																																				
1x11137	Screw Guide 80mm																																																																																				
1x1-1100101	Drill Bit Ø 4.8x180mm Coated - Quick Connect																																																																																				
1x91150	Bone Screw T Wrench																																																																																				
Code	Description																																																																																				
6x93010	Large Clamp																																																																																				
1x932300	Rod Ø 12mm L 300mm																																																																																				
2x932150	Rod Ø 12mm L 150mm																																																																																				
2x912650	Self-drilling XCaliber Screws, L 260mm, thread length 50mm																																																																																				
2x911550	Self-drilling XCaliber Screws, L 150mm, thread length 50mm																																																																																				
1x11138	Drill Guide Ø 4.8mm																																																																																				
1x11137	Screw Guide 80mm																																																																																				
1x1-1100101	Drill Bit Ø 4.8x180mm Coated - Quick Connect																																																																																				
1x91150	Bone Screw T Wrench																																																																																				

9 TI	10 F&A	11 F&A	12 F&A	13 F&A																																																																													
99-93794 Galaxy UNYCO Diaphyseal Tibia Box: <table border="1"> <thead> <tr> <th>Code</th> <th>Description</th> </tr> </thead> <tbody> <tr><td>2x99-93574</td><td>Galaxy UNYCO Mini Kit Tibia Sterile</td></tr> <tr><td>1x99-939250</td><td>Rod Ø 12mm L 350mm Sterile</td></tr> <tr><td>1x99-93509</td><td>Galaxy UNYCO Mini Kit Instruments Sterile</td></tr> </tbody> </table>	Code	Description	2x99-93574	Galaxy UNYCO Mini Kit Tibia Sterile	1x99-939250	Rod Ø 12mm L 350mm Sterile	1x99-93509	Galaxy UNYCO Mini Kit Instruments Sterile	99-93503 Ankle Sterile kit - Transfix Pin Ø6mm 99-93499 Ankle Sterile kit - Transfix Pin Ø4mm: <table border="1"> <thead> <tr> <th>Code</th> <th>Description</th> </tr> </thead> <tbody> <tr><td>4x93010</td><td>Large Clamp</td></tr> <tr><td>2x93030</td><td>Transitional Clamp</td></tr> <tr><td>1x932300</td><td>Rod Ø 12mm L 300mm</td></tr> <tr><td>1x932200</td><td>Rod Ø 12mm L 200mm</td></tr> <tr><td>1x932150</td><td>Rod Ø 12mm L 150mm</td></tr> <tr><td>2x911540</td><td>XCaliber Screws</td></tr> <tr><td>1x11138</td><td>Drill guide d 4.8mm</td></tr> <tr><td>1x11137</td><td>Screw Guide 80mm</td></tr> <tr><td>1x1-1100101</td><td>Drill Bit Ø 4.8x180mm</td></tr> <tr><td>1x91150</td><td>Universal "T" Wrench</td></tr> <tr><td>1x92080</td><td>Transfix Pin thread lenght 60mm, thread Ø 5mm, shaft 4mm</td></tr> </tbody> </table>	Code	Description	4x93010	Large Clamp	2x93030	Transitional Clamp	1x932300	Rod Ø 12mm L 300mm	1x932200	Rod Ø 12mm L 200mm	1x932150	Rod Ø 12mm L 150mm	2x911540	XCaliber Screws	1x11138	Drill guide d 4.8mm	1x11137	Screw Guide 80mm	1x1-1100101	Drill Bit Ø 4.8x180mm	1x91150	Universal "T" Wrench	1x92080	Transfix Pin thread lenght 60mm, thread Ø 5mm, shaft 4mm	99-91647UE XCaliber Ankle Complete Kit, Sterile: <table border="1"> <thead> <tr> <th>Code</th> <th>Description</th> </tr> </thead> <tbody> <tr><td>1x91150</td><td>Universal T Wrench</td></tr> <tr><td>1x91017</td><td>Universal Allen Wrench</td></tr> <tr><td>1x11106</td><td>Drill Guide 3.2x40mm</td></tr> <tr><td>1x11104</td><td>Drill Guide 4.8x40mm</td></tr> <tr><td>1x91047</td><td>XCaliber Ankle Fixator</td></tr> <tr><td>1x11947</td><td>RadioLucent Ankle Pin Guide</td></tr> <tr><td>3x11102</td><td>Screw Guide L 60mm</td></tr> <tr><td>1x91015</td><td>Dynamic Compression/Distraction Unit</td></tr> <tr><td>2x99-611540</td><td>XCaliber Osteotite Screw 150/40</td></tr> <tr><td>3x99-611530</td><td>XCaliber Osteotite Screws 150/30</td></tr> <tr><td>1x1101101</td><td>Cannulated Drill Bit 3.2x200mm</td></tr> <tr><td>2x11014</td><td>X-Wires without olive 1.5x250mm</td></tr> <tr><td>1x1100101</td><td>Drill Bit 4.8x180mm</td></tr> </tbody> </table>	Code	Description	1x91150	Universal T Wrench	1x91017	Universal Allen Wrench	1x11106	Drill Guide 3.2x40mm	1x11104	Drill Guide 4.8x40mm	1x91047	XCaliber Ankle Fixator	1x11947	RadioLucent Ankle Pin Guide	3x11102	Screw Guide L 60mm	1x91015	Dynamic Compression/Distraction Unit	2x99-611540	XCaliber Osteotite Screw 150/40	3x99-611530	XCaliber Osteotite Screws 150/30	1x1101101	Cannulated Drill Bit 3.2x200mm	2x11014	X-Wires without olive 1.5x250mm	1x1100101	Drill Bit 4.8x180mm	99-M1450 Calcaneal Minifixator Kit, sterile: <table border="1"> <thead> <tr> <th>Code</th> <th>Description</th> </tr> </thead> <tbody> <tr><td>1xM145</td><td>Calcaneal External Fixator</td></tr> <tr><td>2x10012</td><td>3mm Allen Wrench</td></tr> <tr><td>1xM210</td><td>T Wrench</td></tr> <tr><td>1xM211</td><td>Driver</td></tr> <tr><td>1xM224</td><td>Calcaneal Template Kit Consisting of: - 6 Wires - 6 Wire Guides - 6 Pin Guides</td></tr> <tr><td>2xM321</td><td>Pin 70/15</td></tr> <tr><td>2xM315</td><td>Pin 70/25</td></tr> <tr><td>4xM317</td><td>Pin 100/30</td></tr> </tbody> </table>	Code	Description	1xM145	Calcaneal External Fixator	2x10012	3mm Allen Wrench	1xM210	T Wrench	1xM211	Driver	1xM224	Calcaneal Template Kit Consisting of: - 6 Wires - 6 Wire Guides - 6 Pin Guides	2xM321	Pin 70/15	2xM315	Pin 70/25	4xM317	Pin 100/30
Code	Description																																																																																
2x99-93574	Galaxy UNYCO Mini Kit Tibia Sterile																																																																																
1x99-939250	Rod Ø 12mm L 350mm Sterile																																																																																
1x99-93509	Galaxy UNYCO Mini Kit Instruments Sterile																																																																																
Code	Description																																																																																
4x93010	Large Clamp																																																																																
2x93030	Transitional Clamp																																																																																
1x932300	Rod Ø 12mm L 300mm																																																																																
1x932200	Rod Ø 12mm L 200mm																																																																																
1x932150	Rod Ø 12mm L 150mm																																																																																
2x911540	XCaliber Screws																																																																																
1x11138	Drill guide d 4.8mm																																																																																
1x11137	Screw Guide 80mm																																																																																
1x1-1100101	Drill Bit Ø 4.8x180mm																																																																																
1x91150	Universal "T" Wrench																																																																																
1x92080	Transfix Pin thread lenght 60mm, thread Ø 5mm, shaft 4mm																																																																																
Code	Description																																																																																
1x91150	Universal T Wrench																																																																																
1x91017	Universal Allen Wrench																																																																																
1x11106	Drill Guide 3.2x40mm																																																																																
1x11104	Drill Guide 4.8x40mm																																																																																
1x91047	XCaliber Ankle Fixator																																																																																
1x11947	RadioLucent Ankle Pin Guide																																																																																
3x11102	Screw Guide L 60mm																																																																																
1x91015	Dynamic Compression/Distraction Unit																																																																																
2x99-611540	XCaliber Osteotite Screw 150/40																																																																																
3x99-611530	XCaliber Osteotite Screws 150/30																																																																																
1x1101101	Cannulated Drill Bit 3.2x200mm																																																																																
2x11014	X-Wires without olive 1.5x250mm																																																																																
1x1100101	Drill Bit 4.8x180mm																																																																																
Code	Description																																																																																
1xM145	Calcaneal External Fixator																																																																																
2x10012	3mm Allen Wrench																																																																																
1xM210	T Wrench																																																																																
1xM211	Driver																																																																																
1xM224	Calcaneal Template Kit Consisting of: - 6 Wires - 6 Wire Guides - 6 Pin Guides																																																																																
2xM321	Pin 70/15																																																																																
2xM315	Pin 70/25																																																																																
4xM317	Pin 100/30																																																																																

14 PAE	15 PAE																												
99-93520 Galaxy Medium - Paediatric Sterile Kit - Ø5mm screw thread: <table border="1"> <thead> <tr> <th>Code</th> <th>Description</th> </tr> </thead> <tbody> <tr><td>2x93140</td><td>Double multipin clamp medium</td></tr> <tr><td>2x939250</td><td>Rod Ø 9mm L250mm</td></tr> <tr><td>4x944540</td><td>Selfdrilling XCaliber Cylindrical Screw shaft Ø 6mm thread 5mm L150/40QC</td></tr> <tr><td>2x11137</td><td>Screw guide 80mm</td></tr> <tr><td>1x30017</td><td>Allen Wrench 5mm</td></tr> <tr><td>1x93160</td><td>QC Wrench</td></tr> </tbody> </table>	Code	Description	2x93140	Double multipin clamp medium	2x939250	Rod Ø 9mm L250mm	4x944540	Selfdrilling XCaliber Cylindrical Screw shaft Ø 6mm thread 5mm L150/40QC	2x11137	Screw guide 80mm	1x30017	Allen Wrench 5mm	1x93160	QC Wrench	99-93521 Galaxy Medium - Paediatric Sterile Kit - Ø4mm screw thread: <table border="1"> <thead> <tr> <th>Code</th> <th>Description</th> </tr> </thead> <tbody> <tr><td>2x93140</td><td>Double multipin clamp medium</td></tr> <tr><td>2x939250</td><td>Rod Ø 9mm L250mm</td></tr> <tr><td>4x944530</td><td>Selfdrilling XCaliber Cylindrical Screw shaft Ø 6mm thread 4mm L150/30QC</td></tr> <tr><td>2x11137</td><td>Screw guide 80mm</td></tr> <tr><td>1x30017</td><td>Allen Wrench 5mm</td></tr> <tr><td>1x93160</td><td>QC Wrench</td></tr> </tbody> </table>	Code	Description	2x93140	Double multipin clamp medium	2x939250	Rod Ø 9mm L250mm	4x944530	Selfdrilling XCaliber Cylindrical Screw shaft Ø 6mm thread 4mm L150/30QC	2x11137	Screw guide 80mm	1x30017	Allen Wrench 5mm	1x93160	QC Wrench
Code	Description																												
2x93140	Double multipin clamp medium																												
2x939250	Rod Ø 9mm L250mm																												
4x944540	Selfdrilling XCaliber Cylindrical Screw shaft Ø 6mm thread 5mm L150/40QC																												
2x11137	Screw guide 80mm																												
1x30017	Allen Wrench 5mm																												
1x93160	QC Wrench																												
Code	Description																												
2x93140	Double multipin clamp medium																												
2x939250	Rod Ø 9mm L250mm																												
4x944530	Selfdrilling XCaliber Cylindrical Screw shaft Ø 6mm thread 4mm L150/30QC																												
2x11137	Screw guide 80mm																												
1x30017	Allen Wrench 5mm																												
1x93160	QC Wrench																												



LEGEND		
● HUMERUS	● ELBOW	● FOREARM
● WRIST		
● PELVIS		
● FEMUR	● KNEE	● TIBIA
● FOOT & ANKLE		
● PAEDIATRICS		



ORDERING INFORMATION

Contact details:

Delivery information:

Manufactured by:
ORTHOFIX Srl
 Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
 Telephone +39 045 6719000, Fax +39 045 6719380

www.orthofix.com 0123

Distributed by:



ORTHOFIX TRAUMA STERILE KITS

HUMERUS | ELBOW | FOREARM

WRIST

PELVIS

FEMUR | KNEE | TIBIA

FOOT & ANKLE

PAEDIATRICS



1 99-93505 Galaxy Shoulder Sterile Kit

2 99-91600UE XCaliber Metadiaphyseal Complete Kit

3 99-93504 Galaxy Elbow Sterile Kit

4 99-93498 Galaxy Small Z-Configuration Sterile Kit

HUMERUS

ELBOW

FOREARM

The modular option with Galaxy Fixation medium components single sterile packed is also

5 99-93601 Galaxy Wrist Sterile Kit

4 99-93498 Galaxy Small Z-Configuration Sterile Kit

6 99-36501 RadioLucent Wrist Fixator Kit, sterile

7 99-93501 Pelvis Sterile Kit

WRIST

PELVIS

8 99-93502 Lower Limb Diaphyseal Sterile Kit

8 99-93502 Lower Limb Diaphyseal Sterile Kit

8 99-93502 Lower Limb Diaphyseal Sterile Kit

9 99-93794 Galaxy UNYCO Diaphyseal Tibia Box

2 99-91600UE XCaliber Metadiaphyseal Complete Kit

FEMUR

KNEE

TIBIA

UNYCO

10 **11** 99-93503 / 99-93499 Ankle Sterile kit

12 99-91647UE XCaliber Ankle Complete Kit, Sterile

13 99-M1450 Calcaneal Minifixator Kit, sterile

FOOT & ANKLE

14 99-93520 Galaxy Medium - Paediatric Sterile Kit - Ø5mm screw thread

15 99-93521 Galaxy Medium - Paediatric Sterile Kit - Ø4mm screw thread

4 99-93498 Galaxy Small Z-Configuration Sterile Kit

PAEDIATRICS

EXTERNAL FIXATION INDICATIONS IN TRAUMA

The use of an external fixator for **damage control** is to provide skeletal stability and for this to be achieved efficiently; that is, **the objective is accomplished quickly and does not hinder future definitive care**. Skeletal stability in damage control places considerable strains on an external fixator system as reduction of the fracture is usually approximate and, as such, bone contact at the fracture site and shared stability with the external fixator may be minimal, thereby **leaving the provision of stability entirely to the external fixator construct**.

The Galaxy Fixation and Galaxy UNYCO have been designed with this characteristic in mind. While Galaxy UNYCO has been specifically developed for damage control and emergency - therefore it is purely for temporary stabilisation - Galaxy Fixation may also be needed for definitive treatment of long bone fractures. The external fixator configuration will depend on the amount of bone contact at the fracture site, the fracture pattern and the segment or segments of bone involved.