



RIVAL™

○ REDUCE

□ VIEW

▶ BITE



VIEW

Plating System

Wide range of
implants and sizes

Optimal implant selection

Designed to

reduce soft tissue irritation

Designed to


reduce screw prominence




ORTHOFIX®

Designed to be used with 2.4mm, 2.7mm and 3.2mm locking and non-locking screws

Locking screws



Non-locking screws



Double lead thread screw design

For fast and controlled insertion




Perimeter loading design

The central window allows:

- monitoring of healing process of fractures
- osteotomies
- joint fusion
- application of additional bone grafts or biologics

Guide wire holes

For easy implant placement during surgery



Low profile





Screws lock flush into plates reducing screw prominence

Universal instrumentation

Designed

- to be used across all RIVAL™ product lines
- to help with inventory management
- to provide a comprehensive range of instruments



Distributed by:



Manufactured by:
 ORTHOFIX Srl
 Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
 Telephone +39 045 6719000, Fax +39 045 6719380



www.orthofix.com

RIVI-1905-PL-E0 A 10/19