




New TL-Rapid Struts Provide Quick Solutions

- Quick and easy frame assembly
- Acute adjustment with micro-compression or distraction
- Lockable universal hinge on both ends of strut



SIMPLE STABLE
VERSATILE

TrueLok is Simple
TrueLok is Stable
TrueLok is Versatile

Short Range	Medium Range	Long Range
70 - 104mm	98 - 188mm	163 - 300mm




Struts are available in three sizes:


Length	Part #
SHORT	50-10170
MEDIUM	50-10180
LONG	50-10190

Bolt allows for quick frame adjustments

Acute adjustment with micro-compression or distraction



Distributed by:



Manufactured by:
ORTHOFIX Srl
Via Delle Nazioni 9, 37012 Bussolengo (Verona), Italy
Telephone +39 045 6719000, Fax +39 045 6719380

CE 0123

www.orthofix.com

TL-1104-PL-E0 D 12/18