





The Power of Modularity


When addressing the needs of patients requiring spinal fusion, options matter. The Firebird Spinal Fixation System is a comprehensive implant and instrument platform designed for degenerative and deformity correction procedures in posterior thoracolumbar fusion applications. The system includes modular and pre-assembled implants and offers multiple head and rod options to address single level to complex multi-level fusions.


Pedicle Screw Options

- Color coded by screw diameter
- Modular and pre-assembled screws
- Optional HA coated modular screws*
- Cannulated and non-cannulated screw options
- Modular bodies
- Top loading body
- Closed body
- Offset body
- Reduction body

Implant Features

- Double lead bone screw thread
- Wide range of screw lengths
- 5.5mm diameter Titanium rods
 - Straight
 - Pre-lordosed






Instruments Designed for Ease of Use

The Firebird Spinal Fixation System is designed with operating room efficiency in mind. The system includes a wide variety of ergonomically designed instruments to meet the individual needs of the surgeon.

- Straight and T-handle options for taps and screw drivers
- 4 unique rod reduction options
- Simple compression and distraction instruments



Interbody Options



Surgeon-Designed Implants

- Bulleted tips for distraction
- Firm grip when inserting & rotating with threaded and rail connection
- Clear implant positioning in situ, with tantulum x-ray markers 1mm from the ends

Ergonomic Designed Instruments

- Firm connection to the implant for an insert and rotate technique
- Low Profile Inserter, ideal for MIS techniques
- Angle choices with the Implant Inserters for ease of insertion
- Monolithic implant trials for efficiency





Straight Spacer

NON-CANNULATED SCREWS

| As | sem | b | ec |
|----|-----|---|----|
| | | | |

| Part# | Description |
|----------------------------|--------------|
| Multi-Axial / Self Tapping | |
| 44-3425 | 4.5mm / 25mm |
| 44-3430 | 4.5mm / 30mm |
| 44-3435 | 4.5mm / 35mm |
| 44-3440 | 4.5mm / 40mm |
| 44-3445 | 4.5mm / 45mm |
| 44-3525 | 5.5mm / 25mm |
| 44-3530 | 5.5mm / 30mm |
| 44-3535 | 5.5mm / 35mm |
| 44-3540 | 5.5mm / 40mm |
| 44-3545 | 5.5mm / 45mm |
| 44-3550 | 5.5mm / 50mm |
| 44-3555 | 5.5mm / 55mm |
| 44-3630 | 6.5mm / 30mm |
| 44-3635 | 6.5mm / 35mm |
| 44-3640 | 6.5mm / 40mm |
| 44-3645 | 6.5mm / 45mm |
| 44-3650 | 6.5mm / 50mm |
| 44-3655 | 7.5mm / 40mm |
| 44-3745 | 7.5mm / 45mm |
| 44-3750 | 7.5mm / 50mm |
| 44-3755 | 7.5mm / 55mm |
| 44-3840 | 8.5mm / 40mm |
| 44-3845 | 8.5mm / 45mm |
| 44-3850 | 8.5mm / 50mm |

Modular

| Part# | Decription |
|---------------------------|--------------|
| Bone Screw / Self Tapping | |
| 44-5425 | 4.5mm / 25mm |
| 44-5430 | 4.5mm / 30mm |
| 44-5435 | 4.5mm / 35mm |
| 44-5440 | 4.5mm / 40mm |
| 44-5445 | 5.5mm / 35mm |
| 44-5540 | 5.5mm / 40mm |
| 44-5545 | 5.5mm / 45mm |
| 44-5550 | 5.5mm / 50mm |
| 44-5555 | 5.5mm / 55mm |
| 44-5635 | 6.5mm / 35mm |
| 44-5640 | 6.5mm / 40mm |
| 44-5645 | 6.5mm / 45mm |
| 44-5650 | 6.5mm / 50mm |
| 44-5655 | 7.5mm / 40mm |
| 44-5745 | 7.5mm / 45mm |
| 44-5750 | 7.5mm / 50mm |
| 44-5755 | 7.5mm / 55mm |
| 44-5840 | 8.5mm / 45mm |

CANNULATED SCREWS

| Assembled | |
|---|--------------|
| Part# | Description |
| Cannulated / Multi-Axial / Self Tapping | |
| 77-7525 | 5.5mm / 25mm |
| 77-7530 | 5.5mm / 30mm |
| 77-7535 | 5.5mm / 35mm |
| 77-7540 | 5.5mm / 40mm |
| 77-7545 | 5.5mm / 45mm |
| 77-7550 | 5.5mm / 50mm |
| 77-7555 | 5.5mm / 55mm |
| 77-7630 | 6.5mm / 30mm |
| 77-7635 | 6.5mm / 35mm |
| 77-7640 | 6.5mm / 40mm |
| 77-7645 | 6.5mm / 45mm |
| 77-7650 | 6.5mm / 50mm |
| 77-7655 | 6.5mm / 55mm |
| 77-7740 | 7.5mm / 40mm |
| 77-7745 | 7.5mm / 45mm |
| 77-7750 | 7.5mm / 50mm |
| 77-7755 | 7.5mm / 55mm |
| 77-7840 | 8.5mm / 40mm |
| 77-7845 | 8.5mm / 45mm |
| 77-7850 | 8.5mm / 50mm |

| Modular | |
|--|--------------|
| Part# | Description |
| Cannulated / Bone Screw / Self Tapping | |
| 77-8535 | 5.5mm / 35mm |
| 77-8540 | 5.5mm / 40mm |
| 77-8545 | 5.5mm / 45mm |
| 77-8550 | 5.5mm / 50mm |
| 77-8555 | 5.5mm / 55mm |
| 77-8635 | 6.5mm / 35mm |
| 77-8640 | 6.5mm / 40mm |
| 77-8645 | 6.5mm / 45mm |
| 77-8650 | 6.5mm / 50mm |
| 77-8655 | 6.5mm / 55mm |
| 77-8740 | 7.5mm / 40mm |
| 77-8745 | 7.5mm / 45mm |
| 77-8750 | 7.5mm / 50mm |
| 77-8755 | 7.5mm / 55mm |
| 77-8840 | 8.5mm / 40mm |
| 77-8845 | 8.5mm / 45mm |

HA COATED SCREWS

| Modular | |
|--|--------------|
| Part# | Decription |
| Cannulated / Bone Screw / Self Tapping / Sterile | |
| 77-5525 | 5.5mm X 25mm |
| 77-5530 | 5.5mm X 30mm |
| 77-5535 | 5.5mm X 35mm |
| 77-5540 | 5.5mm X 40mm |
| 77-5545 | 5.5mm X 45mm |
| 77-5550 | 5.5mm X 50mm |
| 77-5555 | 5.5mm X 55mm |
| 77-5625 | 6.5mm X 25mm |
| 77-5630 | 6.5mm X 30mm |
| 77-5635 | 6.5mm X 35mm |
| 77-5640 | 6.5mm X 40mm |
| 77-5645 | 6.5mm X 45mm |
| 77-5650 | 6.5mm X 50mm |
| 77-5655 | 6.5mm X 55mm |
| 77-5725 | 7.5mm X 25mm |
| 77-5730 | 7.5mm X 30mm |
| 77-5735 | 7.5mm X 35mm |
| 77-5740 | 7.5mm X 40mm |
| 77-5745 | 7.5mm X 45mm |
| 77-5750 | 7.5mm X 50mm |
| 77-5755 | 7.5mm X 55mm |
| 77-5825 | 8.5mm X 25mm |
| 77-5830 | 8.5mm X 30mm |
| 77-5835 | 8.5mm X 35mm |
| 77-5840 | 8.5mm X 40mm |
| 77-5845 | 8.5mm X 45mm |
| 77-5850 | 8.5mm X 50mm |
| 77-5855 | 8.5mm X 55mm |

OTHER IMPLANTS

| Part# | Description |
|----------------|-------------|
| Rod / Straight | |
| 52-2040 | 40mm |
| 52-2050 | 50mm |
| 52-2060 | 60mm |
| 52-2070 | 70mm |
| 52-2080 | 80mm |
| 52-2090 | 90mm |
| 52-2100 | 100mm |
| 52-2110 | 110mm |
| 52-2120 | 120mm |
| 52-2140 | 140mm |
| 52-2160 | 160mm |
| 52-2180 | 180mm |
| 52-2200 | 200mm |
| 52-2450 | 450mm |

Rods, Pre-Lordosed

| 52-6035 | 35mm |
|---------|------|
| 52-6040 | 40mm |
| 52-6045 | 45mm |
| 52-6050 | 50mm |
| 52-6055 | 55mm |
| 52-6060 | 60mm |
| 52-6065 | 65mm |
| 52-6070 | 70mm |
| 52-6075 | 75mm |
| 52-6080 | 80mm |
| | |

| Set Screw | | |
|-----------|--|--|
| 44-2001 | | |
| | | |

Set Screw / Offset Body 44-2002

| Body / Top Loading | | |
|--------------------|--|--|
| 44-2101 | | |
| | | |

Body / Closed

44-2102

| Part# | Description |
|-------------|-------------|
| Offset Body | |
| 44-2308 | 8mm |
| 44-2311 | 11mm |
| 44-2314 | 14mm |

Cross Connectors

Cross Connector / Multi-Axial

| 55-5325 | 25mm |
|---------|------|
| 55-5330 | 30mm |
| 55-5335 | 35mm |
| 55-5340 | 40mm |
| 55-5345 | 45mm |
| 55-5350 | 50mm |
| 55-5355 | 55mm |
| 55-5360 | 60mm |
| 55-5365 | 65mm |
| 55-5370 | 70mm |
| 55-5375 | 75mm |
| 55-5380 | 80mm |
| | |

Cross Connector / Fixed

| 15mm |
|------|
| 17mm |
| 19mm |
| 21mm |
| 23mm |
| 25mm |
| - |

Lateral / Rod to Rod

Lateral Connector / Off Set

| 52-6315 | 15mm |
|---------|------|
| 52-6320 | 20mm |
| 52-6325 | 25mm |

| Rod Connector / Axial | |
|-----------------------|---------------|
| 52-6701 | 5.5mm / 5.5mm |

Rod Connector / Parallel

| 52-6801 | 5.5mm / F-F |
|---------|-------------|
| 52-6805 | 5.5mm / T-T |

REDUCTION IMPLANTS AND INSTRUMENTS

| Part# | Description | |
|---------|---|--|
| 44-2103 | Body, Reduction | |
| 61-0112 | Anti-Splay Cap Remover | |
| 61-0400 | Tab Removal Tool | |
| 61-1265 | Reduction Screw Counter Torque Wrench | |
| 61-1330 | Mono-Axial Reduction Screw Driver | |
| 61-1331 | Multi-Axial Reduction Screw Driver | |
| 61-7530 | 5.5mm x 30mm Mono Axial Reduction Screw, Self-Tapping | |
| 61-7535 | 5.5mm x 35mm Mono Axial Reduction Screw, Self-Tapping | |
| 61-7540 | 5.5mm x 40mm Mono Axial Reduction Screw, Self-Tapping | |

| Part# | Description |
|---------|--|
| 61-7545 | 5.5mm x 45mm Mono Axial Reduction Screw, Self-Tapping |
| 61-7630 | 6.5mm x 30mm Mono Axial Reduction Screw, Self-Tapping |
| 61-7635 | 6.5mm x 35mm Mono Axial Reduction Screw, Self-Tapping |
| 61-7640 | 6.5m m x 40mm Mono Axial Reduction Screw, Self-Tapping |
| 61-7645 | 6.5mm x 45mm Mono Axial Reduction Screw, Self-Tapping |
| 61-7740 | 7.5mm x 40mm Mono Axial Reduction Screw, Self-Tapping |
| 61-7745 | 7.5mm x 45mm Mono Axial Reduction Screw, Self-Tapping |
| 61-7750 | 7.5mm x 50mm Mono Axial Reduction Screw, Self-Tapping |
| 68-0111 | Anti-Splay Cap |

INSTRUMENTS

Case 1, 44-9060

| Part# | Description | Part# | Description |
|---------|------------------------------------|---------|--------------------------|
| 52-1001 | Bone Awl | 52-1026 | 6.5mm Bone Tap, Modular |
| 52-1002 | Bone Probe, Straight, Large | 52-1027 | 7.5mm Bone Tap, Modular |
| 52-1003 | Bone Probe, Curved, Large | 52-1028 | 8.5mm Bone Tap, Modular |
| 52-1004 | Dual Sounder, Small | 52-1830 | Mono-Axial Screw Driver |
| 52-1005 | Dual Sounder, Large | 52-1831 | Multi-Axial Screw Driver |
| 52-1011 | Ratcheting Handle, T-Handle | 52-1832 | Modular Screw Driver |
| 52-1013 | Ratcheting Handle, Straight, Small | 52-1040 | Trial Rod, 90mm |
| 55-1006 | X-ray Marker, Right | 52-1041 | Trial Rod, 200mm |
| 55-1007 | X-ray Marker, Left | 52-1042 | Trial Rod, 450mm |
| 52-1024 | 4.5mm Bone Tap, Modular | 52-1308 | Screw Sizing Template |
| 52-1025 | 5.5mm Bone Tap, Modular | | |

Case 2, 44-9070

| Part# | Description | Part# | Description |
|---------|------------------------------------|---------|---------------------------------------|
| 51-7100 | Hook Holder, Regular, Straight | 52-1061 | Set Screw Driver |
| 52-1013 | Ratcheting Handle, Straight, Small | 52-1101 | Cross-Connector Caliper |
| 52-1038 | Head Adjuster | 52-1102 | Cross-Connector Bender, Right |
| 52-1334 | Decorticating Planar | 52-1103 | Cross-Connector Bender, Left |
| 52-1339 | Multi-Axial Adjustment Driver | 52-1104 | Cross-Connector Driver |
| 52-1581 | Rod Inserter | 52-1512 | Torque Limiting Handle |
| 52-1890 | Hinged Distractor | 52-1765 | Counter Torque Wrench |
| 52-1059 | Set Screw Holder / Driver, Short | | · · · · · · · · · · · · · · · · · · · |
| | | | |

Case 3, 44-9080

| Part# | Description | Part# | Description |
|---------|-----------------------------|---------|---------------------------|
| 51-1989 | Tubular Rod Reducer | 52-1755 | Rod Reducer |
| 51-1990 | Driver, Tubular Rod Reducer | 51-1480 | Rod Gripper |
| 52-1011 | Ratcheting Handle, T-Handle | 52-1070 | In-Situ Rod Bender, Right |
| 52-1046 | Rod Bender | 52-1071 | In-Situ Rod Bender, Left |
| 52-1050 | Rod Pusher | 52-1590 | Distractor |
| 52-1251 | Rocker | 52-1591 | Compressor |







Distributed by:

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician. Proper surgical procedure is the responsibility of the medical professional. Operative techniques are furnished as an informative guideline. Each surgeon must evaluate the appropriateness of a technique based on his or her personal medical credentials and experience. Please refer to the "Instructions for Use" supplied with the product for full information on indications for use, contraindications, warnings, precautions, adverse reactions information and sterilization.



FB-1507 © Orthofix Holdings, Inc. 8/2015

www.orthofix.com

1.888.298.5700