

**LONESTAR**<sup>®</sup>  
CERVICAL STAND ALONE SYSTEM




Serrated Radial Ribs

# LONESTAR®

CERVICAL STAND ALONE SYSTEM






## LONESTAR® Spacer

- Anatomical profile plate, screw and spacer system
- Two axial footprints (15mm x 13mm and 17mm x 15mm)
- Two sagittal profiles (parallel and 7° lordotic)
- 9 heights, 6mm – 14mm (6mm – 9mm standard)
- Tactile and visual confirmation of bone screw locking
- Tantalum marker to determine posterior margin
- Large window for packing graft material
- Enhanced stability design available for 7° lordotic 6mm to 10mm heights




## Variety of Screws

- 3.6mm or 4.0mm diameters
- Self-drilling or self-tapping with square drive
- Lengths 10mm to 18mm (12mm – 16mm standard)

Purple	10mm	
Magenta	12mm	
Gold	14mm	
Dark Blue	16mm	
Green	18mm	

## Screw Angulation

- Screw orientation is 38° cephalad/caudal
- Medial orientation 10° nominal



## Cover Plates

- Tactile and visual confirmation of bone screw locking
- Silver color accommodates 6mm – 9mm implant
- Gold color accommodates 10mm – 14mm implant



CE 0086 **RX Only**

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

**EC REP** Medical Device Safety Service, Schiffgraben 41, 30175, Hannover, Germany

1.888.298.5700  
[www.orthofix.com](http://www.orthofix.com)

LS-1602 © Orthofix Holdings, Inc. 5/2016

 **ORTHOFIX®**